

宁波普瑞供应链服务有限公司

PRIMARY
SCS



**LABORATORY
& MEDICAL**

NINGBO PRIMARY SUPPLY CHAIN SOLUTION CO., LTD.

www.primaryscs.com


LABORATORY

TIPS	01
TEST TUBE	02-03
CENTRIFUGE TUBE & PCR TUBE	03-04
CRYOVIAL TUBE	05-06
RACK & BOX	07-08
PASTEUR PIPETTE	09
CUVETTE CUP	10-11
SPECIMEN CONTAINER	12
EMBEDDING CASSETTE	13
MICROPLATES	14-15
MICROSCOPE SLIDES	16
SLIDE MAILER & BOX	17-18
MICROBIOLOGY	19
PLASTICWARE	20-24
GLASSWARE	25-26
PIPETTES & RACK	27-28
INSTRUMENTS	29-31

MEDICAL

VACUUM TUBE	32
NON-VACUUM TUBE	33
SWAB & CERVICAL SCRAPER & BRUSH	34-35
BLOOD LANCET	36
INJECTION PRODUCTION	37-38
OXYGEN MASK	39-40
DIGITAL THERMOMETER	41
VAGINAL SPECULUM	42
SURGICAL BLADES & SCISSORS	43
STETHOSCOPE	44-46
SPHYGMOMANOMETER	47
BANDAGES	48
MEDICAL ADHESIVE TAPE	49
TEST CASSETTE & URINE BAG	50
GAUZE & COTTON & NON-WOVEN PRODUCTS	51-54
OTHERS	55-57
SHARPS CONTAINER	58-61


Tips



Cod.	Description	Specification	Color
PL0101	Tips for Eppendorf Pipette	5ml	Neutral
PL0102		1000ul with graduation	Blue
PL0103		1000ul	Blue
PL0104		200ul	Yellow
PL0105		300ul	Neutral
PL0106		10ul	Neutral
PL0107	Tips for Eppendorf Pipette with filter	1000ul with graduation	Blue
PL0108		1000ul	Blue
PL0109		200ul	Yellow
PL0110		10ul	Neutral
PL0111	Tips for Gilson Pipette	5ml	Neutral
PL0112		200ul long type	Neutral
PL0113		1250ul	Blue
PL0114		1000ul	Blue
PL0115		200ul	Neutral
PL0116		200ul with graduation	Yellow
PL0117		10ul	Neutral
PL0118	Tips for Gilson Pipette with filter	200ul long type	Neutral
PL0119		1000ul	Blue
PL0120		200ul	Yellow
PL0121		200ul with graduation	Yellow
PL0122		10ul	Neutral

Cod.	Description	Specification	Color
PL0123	Tips for Finland Pipette	5ml	Neutral
PL0124		1250ul	Blue
PL0125		1000ul	Neutral
PL0126		200ul	Neutral
PL0127		200ul with graduation	Neutral
PL0128	Tips for Qiujiang-III	10ul	Neutral
PL0129		4x49mm	Yellow
PL0130	Tips for Qiujiang-IV	6x50mm	Neutral
PL0131	Tips for Qiujiang-III / IV	7.4x70mm	Neutral
PL0132	Tips for DRAGONLAB	10ul	Neutral
PL0133		20ul	Neutral
PL0134		1000ul	Neutral
PL0135	Tips for Kodak Multi-use		Neutral
PL0136	Tips for Finn (Lengthened)	200ul, Lengthened	Yellow
PL0137	Tips for MLA	1000ul	Neutral
PL0138	Tips for Oxford	200ul	Neutral
PL0139	Tips for Mla	200ul	Neutral


Tips with Rack



Cod.	Description	Specification	Material
PL0140	Rack for tips	1000ul / 100wells	PP
PL0141		1000ul / 60wells	PP
PL0142		20ul / 96wells	PP
PL0143		10ul / 96wells	PP
PL0144	Tips with rack(without filter)	1000ul / 100well	PP
PL0145		1000ul / 60well	PP
PL0146		200ul / 96well	PP
PL0147		10ul / 96well	PP

Cod.	Description	Specification	Material
PL0148	Tips with rack(with filter)	1000ul / 100well	PP
PL0149		1000ul / 60well	PP
PL0150		200ul / 96well	PP
PL0151		10ul / 96well	PP

Test Tube



Cod.	Description	Specification	Material
PL0201	Test Tube	12x55mm	PS
PL0202		12x60mm	PS
PL0203		12x75mm	PS
PL0204		13x75mm	PS
PL0205		13x78mm	PS
PL0206		12x100mm	PS
PL0207		13x100mm	PS
PL0208		15x100mm	PS
PL0209		16x100mm	PS
PL0210		16x102mm	PS
PL0211	Test Tube	12x75mm	PP
PL0212		13x75mm	PP
PL0213		12x100mm	PP
PL0214		15x100mm	PP
PL0215	Test Tube	16x102mm	PP
PL0216		12x75mm	Glass
PL0217		13x75mm	Glass
PL0218		12x100mm	Glass
PL0219		13x100mm	Glass
PL0220		15x100mm	Glass
PL0221		16x100mm	Glass
PL0222		16x102mm	Glass

Cod.	Description	Specification	Material
PL0223	Conical bottom Tube	12ml	PS
PL0224		12ml	PP
PL0225		17.5x120mm, 15ml	PP
PL0226		17.5x120mm, 16ml	PS
PL0227	Plastic Tube	15x100mm	PE
PL0228		12x100mm	PE
PL0229		10x100mm	PE
PL0230		12x75mm	PE
PL0231		10x75mm	PE
PL0232	Urine Test Tube	10ml	PE
PL0233	Glass tubes with Aluminium screw cup	28x85mm	Glass
PL0234		18x145mm	Glass
PL0235	Glass tubes with screw cap flat bottom/round bottom	12x75mm	Glass
PL0236		13x100mm	Glass
PL0237		16x100mm	Glass
PL0238		16x150mm	Glass
PL0239		20x150mm	Glass



Cod.	Description	Specification	Material
PL0240	Test Tube (Screw cap, flat bottom)	10ml	PS
PL0241		5ml	PP
PL0242	Test Tube (Round bottom, Screw Cap)	5ml	PP
PL0243	Urine Test Tube	12ml without rim	PS
PL0244		12ml wide mouth	PS
PL0245		12ml screw cap	PS
PL0246		12ml with plug cap	PS
PL0247	Urine Measuring Tube (For Germany KOBOLD urine analyzer)	-	PS
PL0248	Urine measuring tube	10ml	PS
PL0249		10ml	PS
PL0250	Urine collection tube (with transfer pipette)	10ml	PE+PP
PL0251	Extraction Pipe	-	PS/PE

Cod.	Description	Specification	Material
PL0252	HCT tube	Thin	PS
PL0253		Thick	PS
PL0254		12.5cm	PS
PL0255	Tube Cap	Ø12mm	PP
PL0256		Ø13mm	PP
PL0257		Ø16mm	PP
PL0258		Ø12mm	PP
PL0259		Ø13mm	PP
PL0260		Ø16mm	PP
PL0261	Replacing Cap (for Beckman Biochemistry Line)	blue	PE
PL0262	Plug Cap	Ø12mm	PE
PL0263		Ø13mm	PE
PL0264		Ø15mm	PE
PL0265		Ø16mm	PE



Cod.	Description	Specification	Label	Material
Screw Cap Culture Tubes				
PL0266	With Screw Cap	10ml, 16x100mm		PP
PL0267	Unassembled	10ml, 16x100mm	✓	PP

Cod.	Description	Specification	Label	Material
Screw Cap Culture Tubes				
PL0268	With Screw Cap	10ml, 16x100mm		PP
PL0269	Assembled	10ml, 16x100mm	✓	PP

Cod.	Description	Specification	Graduation	Material
Culture Tubes with Dual-position Cap (PE Cap)				
PL0270	Autoclavable (121°C) and	5ml, 12x75mm	With	PP
PL0271	freezable (-80°C), can be centrifuged to 3000g	7ml, 12x100mm	With	PP
PL0272		15ml, 17x100mm	With	PP
PL0273	Autoclavable (121°C) and	5ml, 12x75mm	Without	PP
PL0274	freezable (-80°C), can be centrifuged to 3000g	7ml, 12x100mm	Without	PP
PL0275		15ml, 17x100mm	Without	PP

Cod.	Description	Specification	Graduation	Material
Culture Tubes with Dual-position Cap (PE Cap)				
PL0276	Temperature resistant -20°C to 70°C, can be centrifuged to 1400g	5ml, 12x75mm	With	PS
PL0277		15ml, 17x100mm	With	PS
PL0278		5ml, 12x75mm	Without	PS
PL0279	Temperature resistant -20°C to 70°C, can be centrifuged to 1400g	7ml, 12x100mm	Without	PS
PL0280		14ml, 16x95mm	Without	PS
PL0281		15ml, 17x100mm	Without	PS

Centrifuge Tube


- Made of high transparent, medical PP material with separate or attached cap
- Chemical corrosion and low temperature
- Stable under high centrifuge speed, up to 15000rpm
- Adapted to wide range of temperature from -80°C to 121°C
- E.O, Gamma sterile and non-sterile are available

Micro Centrifuge Tube



Cod.	Bottom	Specification	Description
PL0301	Conical	0.2ml	Snap cap
PL0302	Conical	0.5ml	Snap cap
PL0303	Conical	1.5ml	Snap cap
PL0304	Round	2ml	Snap cap
PL0305	Round	5ml	Snap cap
PL0306	Round	7ml	Snap cap
PL0307	Round	10ml	Snap cap


Cod.	Bottom	Specification	Description
PL0308	Self-Standing	0.5ml	Separate screw cap
PL0309	Self-Standing	1.5ml	Separate screw cap
PL0310	Conical	1.5ml	Separate screw cap
PL0311	Conical	1.5ml	Attached screw cap
PL0312	Conical	0.5ml	Separate screw cap
PL0313	Conical	1.5ml	Snap cap (with color)
PL0314	Round	2.0ml	Snap cap (with color)



Cod.	Bottom	Specification	Graduation	Description
PL0315	Round	10ml	White / Black	Screw cap
PL0316	Conical	10ml	White / Black	Screw cap
PL0317	Round	15ml	White / Black	Screw cap
PL0318	Conical	15ml	White / Black	Screw cap


Cod.	Bottom	Specification	Graduation	Description
PL0319	Self standing	50ml	White / Black	Screw cap
PL0320	Conical	50ml	White / Black	Screw cap
PL0321	Round	50ml	No	Screw cap
PL0322	Round	50ml	No	Snap cap

PCR Tube



Cod.	Description	Vol.	Features	Material	Cap Strips
PL0323	8 Tubes Strips (Detachable)	0.2ml	Thin-Walled, Transparent	PP	Domed
PL0324		0.2ml		PP	Flat
PL0325	12 Tubes Strips (Detachable)	0.2ml	Thin-Walled, Transparent	PP	Domed
PL0326		0.2ml		PP	Flat
PL0327	8 Tubes Strips, Undetachable	0.1ml	Thin-Walled, Transparent	PP	-----
PL0328	8 Tubes Strips, Undetachable	0.1ml	Thin-Walled, White	EU	-----
PL0329	8 Tubes Strips, Detachable	0.1ml	Thin-Walled, White	EU	-----

PCR Plates



Cod.	Description	Wells
PL0330	Semi-skirted PCR Plate	96
PL0331	Unskirted PCR Plate, Raised Well Rim	
PL0332	Unskirted PCR Plate	

Cod.	Description	Wells
PL0333	Skirted PCR Plate	96
PL0334	Unskirted PCR Plate, can be cutted	
PL0335	Unskirted PCR Plate, White	


Cryovial Tube

- Made of Medical grade PP material.
- Maximum RCF: 17000G.
- Cryovial tube with external screw cap is designed for freezing samples, the external screw cap design can reduce the contamination probability during sample treatment.
- Cryovial tube with internal screw cap is designed for freezing samples in gas situation of liquid nitrogen. The silicone gel o-ring can enhance the sealing performance of the tube.
- Large white writing area allows easy marking. 2D Cryovial tube is also supplied, laser-etched 2D barcode on bottom.
- Gamma radiation sterile available.



Cod.	Bottom	Specification	Bottom	Resistant Temperature
PL0401	External Thread	1.5ml Ø10x47mm	Self-Standing	-196°C~-121°C
PL0402	External Thread	2.0ml Ø10x47mm		
PL0403	External Thread	0.5ml Ø10x47mm		
PL0404	Internal Thread	1.0ml Ø12x41mm		
PL0405	External Thread	1.0ml Ø12x41mm		

Cod.	Bottom	Specification	Bottom	Resistant Temperature
PL0406	External Thread	5.0ml	Self-Standing	-196°C~-121°C
PL0407	External Thread	1.8ml		
PL0408	External Thread	1.0ml		
PL0409	Internal Thread	0.2ml	Conical	
PL0410	External Thread	0.2ml		



Cod.	Description	Specification	Bottom	Resistant Temperature
PL0411	External Thread	3.5ml Ø12.5x77mm	Self-Standing	-196°C~-121°C
PL0412	External Thread	3.5ml Ø12.5x75mm	Round	
PL0413	Internal Thread	3.5ml Ø12.5x79mm	Self-Standing	-196°C~-121°C
PL0414	Internal Thread	3.5ml Ø12.5x77mm	Round	
PL0415	Internal Thread	4.5ml Ø12.5x93mm	Self-Standing	
PL0416	Internal Thread	4.5ml Ø12.5x91mm	Round	

2D Cryovial Tube with Concave Cap



Cod.	Description	Specification	Bottom	Resistant Temperature	Color	Packing Info.
2D Cryovial Tube with Concave Cap, TYPE A						
PL0417	2D, Cryo box rack	0.5ml Ø8.5x32mm	Round & Self-Standing	-196°C~-121°C	Orange / Yellow / Green / Blue	100pcs/cryo box, 3200pcs/ctn
PL0418	2D, Bag pack	0.5ml Ø8.5x32mm				500pcs/box, 5000pcs/ctn
PL0419	Non-2D, Cryo box rack	0.5ml Ø8.5x32mm				100pcs/cryo box, 3200pcs/ctn
PL0420	Non-2D, Bag pack	0.5ml Ø8.5x32mm				500pcs/box, 5000pcs/ctn
PL0421	2D, Cryo box rack	1.5ml Ø8.5x47mm				100pcs/cryo box, 2400pcs/ctn
PL0422	2D, Bag pack	1.5ml Ø8.5x47mm				500pcs/box, 5000pcs/ctn
PL0423	Non-2D, Cryo box rack	1.5ml Ø8.5x47mm				100pcs/cryo box, 2400pcs/ctn
PL0424	Non-2D, Bag pack	1.5ml Ø8.5x47mm				500pcs/box, 5000pcs/ctn

2D Cryovial Tube with Flat Cap, TYPE B						
PL0425	2D, Cryo box rack	0.5ml Ø8.5x29mm	Round & Self-Standing	-196°C~-121°C	Orange / Yellow / Green / Blue / White	100pcs/cryo box, 3200pcs/ctn
PL0426	2D, Bag pack	0.5ml Ø8.5x29mm				500pcs/box, 5000pcs/ctn
PL0427	Non-2D, Cryo box rack	0.5ml Ø8.5x29mm				100pcs/cryo box, 3200pcs/ctn
PL0428	Non-2D, Bag pack	0.5ml Ø8.5x29mm				500pcs/box, 5000pcs/ctn
PL0429	2D, Cryo box rack	1.5ml Ø8.5x44mm				100pcs/cryo box, 2400pcs/ctn
PL0430	2D, Bag pack	1.5ml Ø8.5x44mm				500pcs/box, 5000pcs/ctn
PL0431	Non-2D, Cryo box rack	1.5ml Ø8.5x44mm				100pcs/cryo box, 2400pcs/ctn
PL0432	Non-2D, Bag pack	1.5ml Ø8.5x44mm				500pcs/box, 5000pcs/ctn




Cod.	Description	Specification	Bottom	Resistant Temperature
PL0433	Internal Thread	1ml Ø12.3x40mm	Round	-196°C~121°C
PL0434		1ml Ø12.3x41mm	Self-standing, TYPE A	
PL0435		1ml Ø12.3x41mm	Self-standing, TYPE B	
PL0436		2ml Ø12.3x46.5mm	Round	
PL0437		2ml Ø12.3x49mm	Self-standing, TYPE A	
PL0438		2ml Ø12.3x49mm	Self-standing, TYPE B	
PL0439		4ml Ø12.3x69.3mm	Round	
PL0440		4ml Ø12.3x70.3mm	Self-standing, TYPE A	
PL0441		4ml Ø12.3x70.3mm	Self-standing, TYPE B	
PL0442	External Thread	2ml Ø12.4x46.5mm	Round	-196°C~121°C
PL0443		2ml Ø12.4x48.5mm	Self-standing, TYPE A	
PL0444		5ml Ø12.4x89.2mm	Round	
PL0445		5ml Ø12.4x90.4mm	Self-standing, TYPE A	

Cod.	Description	Specification	Bottom	Resistant Temperature	Sterilization
PL0446	External Thread without graduation	0.5ml Ø11.2x30mm	Round	-80°C~121°C	No
PL0447		0.5ml Ø11.2x31mm	Self-Standing		No
		0.5ml Ø11.2x31mm	Self-Standing		Gamma
PL0448		1.5ml Ø11.2x47.7mm	Self-Standing		No
		1.5ml Ø11.2x47.7mm	Self-Standing		Gamma
	External Thread with graduation	1.5ml Ø11.2x47.7mm	Self-Standing		Gamma
		1.5ml Ø11.2x47.7mm	Self-Standing		No

Rack for Centrifuge Tube


- Adapted to wide range temperature from -70°C to 121°C
- Reusable after clean or autoclave sterile
- Multicolor can be customized
- Non-sterile



Cod.	Description	Specification	Material
PL0501	Rack for Centrifuge Tube	0.2ml/96wells	PP
PL0502		0.5ml/72wells	PP
PL0503		1.5ml/50wells	PP
PL0504		1.5ml/72wells	PP
PL0505		1.5ml/100wells	PP
PL0506		0.5ml/50wells	PP
PL0507		1.5ml/50wells	PP

Cod.	Description	Specification	Material
PL0508	Rack for Centrifuge Tube	0.5ml, 1.5ml/60wells	PP
PL0509		0.5ml, 1.5ml/96wells	PP
PL0510		15ml/42wells	Cardboard
PL0511		50ml/16wells	Cardboard
PL0512	Multi-Purpose Rack	50ml/21wells	PP
PL0513		15ml/50ml	PP
PL0514		15ml/50ml	PP

Box for Cryovial Tube



Cod.	Description	Specification	Material
PL0515-1	Box for Cryovial Tube	1.8ml / 50wells	PP
PL0515-2		1.8ml / 81wells	PP
PL0515-3		1.8ml / 100wells	PP
PL0516-1	Box for Cryovial Tube	1.8ml / 81wells	PP
PL0516-2		1.8ml / 100wells	PP
PL0517	Box for Cryovial Tube attached cap	1.8ml / 100wells	PP


Cod.	Description	Specification	Material
PL0518-1	Box for Cryovial Tube individual cap	1.8ml / 81wells	PP
PL0518-2		1.8ml / 100wells	PP
PL0519-1	Box for Cryovial Tube	1.8ml / 81wells	Cardboard
PL0519-2		1.8ml / 100wells	Cardboard
PL0519-3		5ml / 36wells	Cardboard

Other Racks



Cod.	Description	Specification
PL0520	Rack for Serological Pipette	Ø13mm
PL0521	Blood sedimentation rack	10-position

Cod.	Description	Specification
PL0522	ESR fast detector	10-position
PL0523	ESR pipette rack	10 pcs



Cod.	Description	Specification	Material
PL0524-1	Rack for Test Tube	Ø13mm / 50wells	PS+PP
PL0524-2		Ø15mm / 50wells	PS+PP
PL0524-3		Ø18mm / 40wells	PS+PP
PL0525-1	Rack for Test Tube	Ø16mm / 50wells	ABS
PL0525-2		Ø18mm / 50wells	ABS
PL0526-1	Rack for Test Tube	Ø13mm	PP
PL0526-2		Ø15mm	PP
PL0526-3		Ø18mm	PP
PL0527-1	Multi-Purpose Test Tube Rack	25wells	ABS
PL0527-2		50wells	ABS

Cod.	Description	Specification	Material
PL0528	Multi-Purpose Test Tube Rack, Orange	50wells	PP
PL0529-1	Rack for Test Tube Z Type	Ø13mm / 25wells	Stainless Steel
PL0529-2		Ø13mm / 50wells	Stainless Steel
PL0529-3		Ø17mm / 25wells	Stainless Steel
PL0529-4		Ø17mm / 50wells	Stainless Steel
PL0530-1	Rack for Test Tube	Ø13mm	Stainless Steel
PL0530-2		Ø15mm	Stainless Steel
PL0530-3		Ø18mm	Stainless Steel




Cod.	Description	Specification
PL0531	Glass Pasteur Pipette	230mm
PL0532		150mm

Serological Pipette



Cod.	Spec.	Negative Graduation	Color
PL0533	1ml	-0.4ml	Yellow
PL0534	2ml	-0.4ml	Green
PL0535	5ml	-3.0ml	Blue
PL0536	10ml	-3.0ml	Orange
PL0537	25ml	-3.0ml	Red
PL0538	50ml	-5.0ml	Black

ESR Pipette System (Westergren)



Cod.	Spec.
PL0539-1	ESR Pipette
PL0539-2	ESR Pipette (with 3.8% sodium citrate diluents)
PL0540-1	ESR Pipette
PL0540-2	ESR Pipette (with 3.9% sodium citrate diluents)

Cod.	Spec.
PL0541	ESR Pipette
PL0542	ESR Pipette
PL0543	ESR Pipette
PL0544	ESR Rack, 10-Position

Pasteur Pipette

- Made of superior LDPE material, Non-toxic, excellent transparency
- Various sizes to choose, with molded graduation or not
- Available in bulk or individual pack
- E.O. Gamma sterile and non-sterile are available



Cod.	Description	Length	Overall Vol.	Drop Vol.
PL0601	Pasteur Pipette	65mm	--	10ul
PL0602		70mm	--	15ul
PL0603		80mm	--	20ul
PL0604		148mm	--	15ul
PL0605		67mm	--	20ul
PL0606		65mm	--	48ul
PL0607		66mm	0.2ml	48ul
PL0608		85mm	--	40ul
PL0609		115mm	--	46ul
PL0610		130mm	0.5ml	40ul
PL0611		160mm	--	40ul
PL0612		145mm	--	45ul
PL0613		230mm	3ml	40ul
PL0614		350mm	10ml	55ul
PL0615		143mm	1ml	42ul
PL0616		145mm	1ml	45ul
PL0617		160mm	1ml	42ul
PL0618		150mm	2ml	48ul
PL0619		160mm	3ml	48ul
PL0620		162mm	3ml	52ul

Cod.	Description	Length	Overall Vol.	Drop Vol.
PL0621	Double-bulb Transfer Pipette	64mm	80ul	25ul
PL0622		85mm	100ul	30ul
PL0623	Pasteur Pipette	78mm	--	35ul
PL0624		110mm	--	40ul
PL0625		182mm	3ml	48ul
PL0626		132mm	--	48ul
PL0627		122mm	--	30ul
PL0628		38mm	--	5ul
PL0629		75mm	--	50ul
PL0630		65mm	--	20ul
PL0631		80mm	--	20ul
PL0632		101mm	--	30ul
PL0633		103mm	--	30ul
PL0634	Double-bulb Transfer Pipette	115mm	--	40ul
PL0635	Pasteur Pipette	118mm	--	40ul
PL0636		130mm	--	30ul



Cod.	Description	Material
PL0701	For Spectrophotometers, 4.5ml Micro	PS
PL0702	For Spectrophotometers, 2.5ml Semi-micro	PS
PL0703	For Spectrophotometers, 2.5ml Semi-micro	PS
PL0704	For German BE Coagulometer 1 Channel	PS
PL0705	For German BE Coagulometer 2 Channel	PS
PL0706	For German BE Coagulometer 4 Channel	PS
PL0707	For German BE Automatic	Black Frame ABS + PP
PL0708	For German BE Automatic 8 Channel	Black Frame ABS + PP
PL0709	For German BECL4 2 Channel	PS
PL0710	For German BECL4 4 Channel	PS
PL0711	For German BECL4 8 Channel	PS
PL0712	For German TECO	PS
PL0713	For German TECO 2 Channel	PS
PL0714	For German TECO 1 Channel, size S	PS
PL0715	For German TECO 1 Channel, size L	PS

Cod.	Description	Material
PL0716	For Sysmex CA 530	PS
PL0717	For Sysmex CA 1500	PS
PL0718	For Sysmex CA 6000	PS
PL0719	For Sysmex CA 7000	PS
PL0720	For Spain Cruor Apparatus	PS
PL0721	For German LABOR	PS
PL0722	For American Pacific, German TECO 4 Channel	PS
PL0723	For Beckman 8 Channel	PS
PL0724	For Precil & French STAGO	PS
PL0725	For Steellex Coagulometer	PS
PL0726	For American Coulter Coagulometer	PS
PL0727	For Beijing Perlong Coagulometer	PS
PL0728	For Swiss SEROZYME Coagulometer	PS
PL0729	For Roche Coagulometer	PS

Cuvette Cup(Sample Cup)



Cod.	Description	Material
PL0730	For Roche Immunoassay Analyzer	PS
PL0731	For Roche Immunoassay Analyzer (Rack Packed)	PS
PL0732	Tips for Roche Immunoassay Analyzer	ABS
PL0733	Tips for Roche Immunoassay Analyzer (Rack Packed)	ABS
PL0734	For Roche E170, E601 Biochemical Analyzer	PS + ABS
PL0735	For Beckman ACL-TOP	PS


Cod.	Description	Material
PL0736	For Siemens Luminescence Immunoassay Analyzer	PP
PL0737	For Beckman Immunoassay System 81901	PP
PL0738	For LABiTec CoaLAB1000 Coagulometer	PS
PL0739	For French STAGO Fully Auto Coagulometer	PS
PL0740	For Hitachi 902 / 7020 Coagulometer	PS + ABS
PL0741	For Hitachi 904 / 7060 Coagulometer	PS + ABS



Cod.	Description	Specification	Material
PL0742	For Swiss AC900 series	Ø24x45	PS + PP
PL0743	For Sysmex series	Ø30x70	PE
PL0744	For Japan Kohden / Botest series	Ø30x60	PE
PL0745	For Mindray BC-2000	Ø30x52	PE
PL0746	For Hitachi 7150, 7060 series	Ø16x38	PS
PL0747	For American Beckman CX series	Ø13x24	PS
PL0748	For American Beckman 700 & American Dupont	Ø14x25	PS
PL0749	For Kodak Biochemical Analyzer	Ø13x20	PS
PL0750	For American IL	Ø13x17	PS
PL0751	For American Corning 550, 560 series	Ø12x20	PS
PL0752	For Roche COBAS	Ø8x34	PP
PL0753	For Olympus AU	Ø13x38	PP

Cod.	Description	Specification	Material
PL0754	For American Abbott CCX	Ø15x40	PS
PL0755	For Hitachi Minimum Analyzer	Ø13x38	PS
PL0756	For Bowlinman / Roche	Ø7x26	PS
PL0757	For Shanghai Antai	Ø10x19	PP
PL0758	For American M10	Ø11x19	PP
PL0759	For SINNOWA	--	PP
PL0760	For Hitachi (Conical sample cup)	--	PS
PL0761	For American Abbott	--	PS
PL0762	For American Abbott AXSYM	Ø17.5x53	PS
PL0763	For Ciba Corning 550 Express	--	PS
PL0764	For Finland Cornyi	--	PS
PL0765	For American MD2000 / 1800	--	PS
PL0766	For American MD3000	--	PS

Cuvette Cup(Sample Cup)



Cod.	Description	Material
PL0767	For Abbott Architect i2000 SR	PP
PL0768	For Simens	PP
PL0769	For MEDICA Ø13x24	PS
PL0770	For Beckman DXI800	PP
PL0771	For DiaSorin biochemical Cup (six holes)	PS
PL0772	For Dade Behring protein analyzer	PS
PL0773	Quantitative Suction Dropper, 5ul	Dropper: PS Tube & Cap: PP
PL0774	Quantitative Suction Dropper, 10ul	
PL0775	Quantitative Suction Dropper, 7.5ul	

Cod.	Description	Material
PL0776	For Rayto Coagulometer	PS
PL0777	For SEKISUI COAPRESTA 2000	PS
PL0778	For German Sigma KC-40, size S	PS
PL0779	For Thermo Fisher Indiko & Plus Analyzer	PS
PL0780	For MOL-200L Biochemical Analyzer	PS
PL0781	For Abbott ALCYON Biochemical Analyzer	PS



Cod.	Vol.	Cap	Cup Material
PL0801	120ml	Screw cap	PP
PL0802	120ml	Screw cap	PP
PL0803	120ml	Screw cap, concave cover	PP
PL0804	100ml	Screw cap	PP
PL0805	90ml	Screw cap, concave cover	PP
PL0806	60ml	Screw cap	PP
PL0807	60ml	Screw cap	PP
PL0808	40ml	Screw cap	PS
PL0809	40ml	Screw cap	PP

Cod.	Vol.	Cap	Cup Material
PL0810	40ml	Screw cap	PP
PL0811	40ml	Snap cap + Label	PP
PL0812	40ml	---	PVC (Transparent)
PL0813	40ml	---	PVC (White)
PL0814	30ml	Screw cap	PS / PP
PL0815	30ml	Screw cap with Leak-Proof mat	PS
PL0816	30ml	Screw cap with Leak-Proof mat	PP
PL0817	30ml	Snap cap, with cover	PP




Cod.	Vol.	Cap	Cup Material
PL0818	30ml	Snap cap	PP
PL0819	25ml	Screw cap	PS
PL0820	25ml	Screw cap	PP
PL0821	20ml	Screw cap	PS
PL0822	20ml	Screw cap	PS
PL0823	20ml	Snap cap	PP

Stool Container




Cod.	Vol.	Cap	Cup Material
PL0824	60ml	Screw cap	PP
PL0825	60ml	Screw cap	PP
PL0826	60ml	Screw cap	PP
PL0827	40ml	Screw cap	PS
PL0828	40ml	Screw cap	PP
PL0829	40ml	Snap cap	PP (White)
PL0830	40ml	Snap cap	PP
PL0831	30ml	Snap cap	PP (Black)
PL0832	30ml	Snap cap	PP (White)
PL0833	30ml	Screw cap	PP

Cod.	Vol.	Cap	Cup Material
PL0834	30ml	Screw cap	PP / PS
PL0835	30ml	Screw cap	PP
PL0836	25ml	Snap cap	PP
PL0837	10ml	Snap cap	PS
PL0838	10ml	Snap cap	PS
PL0839	5ml	Snap cap	PP
PL0840	5ml	Snap cap	PP
PL0841	In anus type	Snap cap	PS
PL0842	Stool Stick	--	PE




Cod.	Description	Holes size	Writing angle	Material
PL0901	One-piece & snap-latch	1.6mm x 1.6mm	45°	POM
PL0902	One-piece & snap-latch	0.9mm x 0.9mm	45°	POM
PL0903	Special hinge, Unassembled	0.9mm x 0.9mm	45°	POM
PL0904	Special hinge, Assembled	0.9mm x 0.9mm	45°	POM
PL0905	Two-pieces, Unassembled	0.9mm x 0.9mm	45°	POM
PL0906	Two-pieces, Assembled	0.9mm x 0.9mm	45°	POM
PL0907	Two-pieces, one compartment, Assembled	0.26mm x 0.26mm	45°	POM
PL0908	Two-pieces, four compartments, Assembled	0.26mm x 0.26mm	45°	POM
PL0909	One-piece & snap-latch	5.0mm x 1.0mm	45°	POM
PL0910	One-piece & snap-latch	5.0mm x 1.0mm	35°	POM
PL0911	One-piece, with extra deflector holes	5.0mm x 1.0mm	45°	POM
PL0912	Special hinge, Unassembled	5.0mm x 1.0mm	45°	POM
PL0913	Special hinge, Assembled	5.0mm x 1.0mm	45°	POM
PL0914	Two-pieces, one compartment, Unassembled	5.0mm x 1.0mm	45°	POM
PL0915	Two-pieces, one compartment, Assembled	5.0mm x 1.0mm	45°	POM
PL0916	Two-pieces, four compartments, Unassembled	5.0mm x 1.0mm	45°	POM
PL0917	Two-pieces, four compartments, Assembled	5.0mm x 1.0mm	30°	POM
PL0918	One-piece	5.0mm x 1.0mm	30°	POM
PL0919	Two-pieces, Unassembled	5.0mm x 1.0mm	30°	POM
PL0920	Two-pieces, Assembled	5.0mm x 1.0mm	30°	POM
PL0921	Cassette, without lid	Ø2.4mm	45°	POM
PL0922	Cassette, without lid	5.0mm x 1.0mm	45°	POM
PL0923	Cassette, without lid	1.6mm x 1.6mm	45°	POM



Cod.	Description	Specification	Material
PL0924	Stainless Steel Lids	45°	Stainless steel
PL0925	Embedding Rings	45°	ABS
PL0926	Biopsy Foam Pads	30x25x2mm	PU Foam

Cod.	Description	Net Well Dim.	External Dim.	Material
PL0927-1	Stainless steel base molds	7x7x7mm	52x35x11mm	Stainless steel
PL0927-2		15x15x7mm	52x35x11mm	Stainless steel
PL0927-3		24x24x7mm	52x35x11mm	Stainless steel

Cod.	Description	Net Well Dim.	External Dim.	Material
PL0927-4	Stainless steel base molds	30x24x7mm	52x35x11mm	Stainless steel
PL0927-5		37x24x7mm	52x35x11mm	Stainless steel
PL0928-1	PVC base molds	7x7x5mm	50x37x12mm	PVC
PL0928-2		17x17x5mm	50x37x12mm	PVC
PL0928-3		25x25x5mm	50x37x12mm	PVC
PL0928-4		25x32x5mm	50x37x12mm	PVC
PL0928-5		25x40x5mm	50x37x12mm	PVC
PL0928-6		30x20x12mm	55x38x18mm	PVC




Cod.	Description	Specification	Material	Treatment
PL1001	Tissue Culture Plate	6wells	PS	TC treated
PL1002		12wells	PS	
PL1003		24wells	PS	
PL1004		48wells	PS	
PL1005		96wells, U-bottom	PS	
PL1006		96wells, V-bottom	PS	
PL1007		96wells, flat-bottom	PS	
PL1008		96wells, flat-bottom	PS (black)	
PL1009		96wells, flat-bottom	PS (white)	
PL1001-1	Tissue Culture Plate	6wells	PS	Non-treated
PL1002-1		12wells	PS	
PL1003-1		24wells	PS	
PL1004-1		48wells	PS	
PL1005-1		96wells, U-bottom	PS	
PL1006-1		96wells, V-bottom	PS	
PL1007-1		96wells, flat-bottom	PS	
PL1008-1		96wells, flat-bottom	PS (black)	
PL1009-1		96wells, flat-bottom	PS (white)	

Cod.	Description	Specification	Material	Treatment
PL1010-1	Demountable Tissue Culture Plate	96wells, flat-bottom	HIPS+PS	--
PL1010-2		96wells, U-bottom	HIPS+PS	
PL1010-3		48wells, U-bottom	HIPS+PS	
PL1011	96wells Detachable Elisa Plate	12 wells detachable strip	PS	--
PL1012		8 wells detachable strip	PS	
PL1013		8 wells detachable strip	PS	
PL1014		8wells detachable strip +Concave sealing cap	PS	



Cod.	Description	Specification	Diameter / Depth	Material
PL1015	96wells Deepwell Plate with Silica Gel Pad	Square well, U-bottom	8.0 / 38.0 mm	PP
PL1016		Square well, U-bottom	8.0 / 20.8 mm	PP
PL1017		Square well, V-bottom	7.2 / 27.2 mm	PP
PL1018		Round well, Flat bottom, Liquid Release Hole	8.3 / 29.0 mm	PP
PL1015-1	96wells Deepwell Plate without Silica Gel Pad	Square well, U-bottom	8.0 / 38.0 mm	PP
PL1016-1		Square well, U-bottom	8.0 / 20.8 mm	PP
PL1017-1		Square well, V-bottom	7.2 / 27.2 mm	PP
PL1018-1		Round well, Flat bottom, Liquid Release Hole	8.3 / 29.0 mm	PP



Cod.	Description	Specification	Material
PL1019-1	Staining Plate	12wells, Ø15-16	Glass
PL1019-2		8wells, Ø15	Glass
PL1019-3		6wells, Ø20-22	Glass
PL1019-4		12wells, Ø20-22	Glass
PL1019-5		24wells, Ø20-22	Glass
PL1019-6		30wells, Ø20-22	Glass

Cod.	Description	Specification	Material
PL1019-7	Staining Plate	24wells, Ø24	Glass
PL1019-8		30wells, Ø24	Glass
PL1020	Hemagglutination Plate	96holes, disposable	PVC
PL1021		96holes	Glass
PL1022	Reagent Groove	V-bottom	PVC

Petri Dish

- Made of high transparency medical grade PS material, used for fungus, bacteria and other microorganism culture
- Easy to stack up with the stack up circle
- Vent design for easy air exchange
- E.O and Gamma sterile are available



Cod.	Specification	Treatment	Material
PL1101	Ø35mm	Non-TC treated	PS
PL1102	Ø60mm		PS
PL1103	Ø60mm, Ring gear		PS
PL1104	Ø70mm		PS
PL1105	Ø90mm		PS
PL1106	Ø90mm, 2rooms		PS
PL1107	Ø90mm, 3rooms		PS
PL1108	Ø90mm, Ring gear		PS
PL1109	Ø90mm, Radian Design		PS
PL1110	100mm Square		PS
PL1111	120mm Square		PS
PL1112	Ø150mm		PS
PL1113	100mm Square with Grid		PS
PL1114	Ø65mm with Grid		PS
PL1115	Ø55mm contact dish with Grid		PS

Cod.	Specification	Treatment	Material
PL1101-1	Ø35mm	TC treated	PS
PL1102-1	Ø60mm		PS
PL1103-1	Ø60mm, Ring gear		PS
PL1104-1	Ø70mm		PS
PL1105-1	Ø90mm		PS
PL1106-1	Ø90mm, 2rooms		PS
PL1107-1	Ø90mm, 3rooms		PS
PL1108-1	Ø90mm, Ring gear		PS
PL1109-1	Ø90mm, Radian Design		PS
PL1110-1	100mm Square		PS
PL1111-1	120mm Square		PS
PL1112-1	Ø150mm		PS
PL1113-1	100mm Square with Grid		PS
PL1114-1	Ø65mm with Grid		PS
PL1115-1	Ø55mm contact dish with Grid		PS



Individual pack, Sterilization



Bulk Pack, Sterilization




Cod.	Description	Volume	Cap style	Treatment	Material	Sterilization
PL1116	Cell Culture Flask	50ml	Standard	TC treated	PS	Gamma Radiation
			Vent		PS	
PL1117		250ml	Standard		PS	
			Vent		PS	
PL1118		600ml	Standard		PS	
			Vent		PS	
PL1116-1		50ml	Standard	Non-TC treated	PS	Gamma Radiation
			Vent		PS	
PL1117-1		250ml	Standard		PS	
			Vent		PS	
PL1118-1		600ml	Standard		PS	
			Vent		PS	

Microscope Slides & Cover Glass

• Designed for routine experiments where no extra adhesion or no adaptability to auto-writers is required in histology-pathology, hematology, cytology, microbiology and etc, for blood or body liquid tests, frozen tissue section examinations, biological specimen preparations and etc..

Microscope Slides



Cod.	Specification
7101	Ground Edges
7101B	Plain, Beveled Edges, 45°
7102	Cut Edges
7103	Concave, Ground Edges
7104	Double concave, Ground Edges
7105	Ground Edges Single frosted end
7105B	Single frosted end Beveled Edges, 45°

Cod.	Specification
7105-1	Cut Edges Single frosted end
7107	Ground Edges Single frosted end
7107B	Twin frosted ends Beveled Edges, 45°
7107-1	Cut Edges Single frosted end
7109	blue/green/pink/orange/yellow/white 90° corners
7109-1	blue/green/pink/orange/yellow/white 45° corners


Cover Glass




Cod.	Specification
7201	10x10mm
	12x12mm
	15x15mm
	18x18mm
	20x20mm
	22x22mm
	24x24mm
	24x36mm
	24x40mm
	24x43mm
	24x50mm
	24x60mm

Cod.	Specification
7201-1	Ø8mm
	Ø9mm
	Ø10mm
	Ø11mm
	Ø12mm
	Ø13mm
	Ø14mm
	Ø16mm
	Ø18mm
	Ø20mm
	Ø22mm

Others




Cod.	Description	Specification
PL1201	Positive Charged Slides	These Color-Frosted(+)glass slides are made by a new process which offers the following features: It places a permanent positive charge in the microscope slide.
PL1202	Counting Chamber	Improved Neubauer Thoma, new Fuchs-Rosenthal
PL1203	Cover Slip for Chamber	22x26mm, thickness 0.5mm




Cod.	Description	Specification	Material
PL1301	Slide Mailer	1 place	PP
PL1302		2 place	PP
PL1303	Slide Mailer	3 place	PP
PL1304		2 place	PP
PL1305	Slide Mailer	2 place	PP
PL1306		3 place	PP
PL1307	Slide Mailer	5 place	PP
PL1308		5 place	PP
PL1309	Slide Storage Box with Cork Pad	5 place	PP
PL1310		25 place	ABS
PL1311		50 places	ABS
PL1312		100 places	ABS

Cod.	Description	Specification	Material
PL1313	Slide Storage Box with Paper Pad	25 place	ABS
PL1314		50 places	ABS
PL1315		100 places	ABS
PL1316	Slide Storage Box with Paper Pad	25 places	Wooden
PL1317		50 places	Wooden
PL1318		100 places	Wooden
PL1319	Slide Dispenser	--	ABS
PL1320	Slide Tray	20 places	ABS



Cod.	Description	Specification	Material
PL1318	Staining Jar and Rack	24 places	PP
PL1319	Staining Rack	30 places	Stainless Steel
PL1320		60 places	Stainless Steel
PL1321		20 places	Stainless Steel
PL1322		50 places	PPS
PL1323		30 places	PPS

Cod.	Description	Specification	Material
PL1324	Staining Jar with Lid	5 places	Glass
PL1325		9 places	Glass
PL1326		10 places	Glass
PL1327		20 places	Glass
PL1328		30 places	Glass
PL1329		60 places	Glass



Cod.	Description	Specification	Features	Material
PL1330	Slide Mailer	1 place	with lid and dividers	Cardboard
PL1331		2 place	with lid and dividers	
PL1332		3 place	with lid and dividers	
PL1333		4 place	with lid and dividers	
PL1334		9 place	without lid or dividers	
PL1335		10 place	with lid and dividers	
PL1336		20 place	with lid and dividers	



Cod.	Description	Specification	Features	Material
PL1337	Slide Mailer	20 place	without lid, with dividers	Cardboard
PL1338		24 place	with lid, without dividers	
PL1339		32 place	with lid, without dividers	
PL1340		30 place	with lid and dividers	
PL1341		48 place	with lid and dividers	

Inoculating Loop

- Four types of packaging are available.
- Type A: Individual pack, 5000pcs/ctn

- Type B: 5pcs/pack, 5000pcs/ctn
- Type C: 10pcs/pack, 10000pcs/ctn

- Type D: 20pcs/pack, 10000pcs/ctn
- Available in E.O. or Gamma radiation sterile.




Cod.	Description	Specification	Material	Color	Packing Info.
PL1401	Rigid Loop	Needle	PS	Neutral	A / B / C / D
PL1402		1ul	PS	Blue	A / B / C / D
PL1403		10ul	PS	Violet	A / B / C / D
PL1404	Flexible Loop	Needle	HIPS	White	A / B / C / D
PL1405		1ul	HIPS	Blue	A / B / C / D
PL1406		10ul	HIPS	Yellow	A / B / C / D
PL1407	Rigid Loop	1ul	AS	Neutral	A / B / C / D
PL1408		10ul	AS	Blue	A / B / C / D
PL1409		1ul + 10ul	AS	Yellow	A / B / C / D
PL1410	Flexible Loop	1ul	ABS	White	A / B / C / D
PL1411		10ul	ABS	Blue	A / B / C / D
PL1412		1ul + 10ul	ABS	Yellow	A / B / C / D
PL1413	Nichrome Wire Loop	Length 200mm	--	Black	1pcs/pk, 120pcs/ctn
PL1414		Replacement loop end, $\phi 2.5 \times 50\text{mm}$	Nichrome	--	50pcs/bottle
PL1415		Replacement needle end, length 50mm	Nichrome	--	50pcs/bottle
PL1416		Replacement lancet end, length 50mm	Nichrome	--	50pcs/bottle

Cell Spreader & Scraper & Lifter

- Four types of packaging are available.
- Type A: Individual pack, 3000pcs/ctn

- Type B: 5pcs/pack, 3000pcs/ctn
- Type C: 10pcs/pack, 3000pcs/ctn

- Type D: Bulk pack, 6000pcs/ctn
- Available in E.O. or Gamma radiation sterile.



Cod.	Description	Specification	Material
	Bacterial Cell Spreader		
PL1417	L shape	L:148mm x W:30mm	PP
PL1418-1	Triangle shape	L:215mm x W:30mm	PP
PL1418-2		L:240mm x W:60mm	PP
PL1418-3		L:160mm x W:40mm	ABS
PL1419-1	T shape	L:140mm x W:35mm	PP
PL1419-2		L:182mm x W:32mm	ABS

Cod.	Description	Specification	Material
	Cell Scraper & Lifter		
PL1420-1	Cell Scraper	L:180mm x W:20mm	TPE+ABS
PL1420-2	Cell Scraper	L:180mm x W:30mm	TPE+ABS
PL1421-1	Cell Scraper	L:250mm x W:20mm	TPE+ABS
PL1421-2	Cell Scraper	L:250mm x W:30mm	TPE+ABS
PL1422	Cell Lifter	L:150mm x W:20mm	PE
PL1423-1	Cell Lifter	L:215mm x W:10mm	PE
PL1423-2	Cell Lifter	L:215mm x W:20mm	PE
PL1424-1	Cell Scraper	L:180mm x W:20mm	TPE+ABS
PL1424-2	Cell Scraper	L:180mm x W:30mm	TPE+ABS
PL1425-1	Cell Scraper	L:250mm x W:20mm	TPE+ABS
PL1425-2	Cell Scraper	L:250mm x W:30mm	TPE+ABS

Plastic Beaker

- Suitable for general laboratory use, such as obtaining quantities of solution or containing buffer and work reagent solution for short time
- Transparent, chemical resistance, durable and autoclavable
- No drip, single-spout design provides precision pouring control



Cod.	Specification	Capacity	Material
PL1501	With graduation	25ml	PP
PL1502		50ml	PP
PL1503		100ml	PP
PL1504		150ml	PP
PL1505		250ml	PP
PL1506		300ml	PP
PL1507		500ml	PP
PL1508		1000ml	PP
PL1509	With handle, With graduation	250ml	PP
PL1510		500ml	PP
PL1511		1000ml	PP
PL1512		2000ml	PP
PL1513		3000ml	PP
PL1514		5000ml	PP



Measuring Cylinder

- The base keeps cylinders stable on the bench top preventing tipping hazards-won't break or chip like glassware
- Molded graduation or printed graduation, chemical resistance and break resistance

Cod.	Capacity	Capacity	Material
PL1515	With graduation	50ml	PP
PL1516		100ml	PP
PL1517		250ml	PP
PL1518		500ml	PP
PL1519		1000ml	PP

Bottles and Funnel



Cod.	Specification	Capacity	Material
PL1520-1	Conical Flask with cap	250ml	PP
PL1520-2		500ml	PP
PL1521-1	Washing bottle with red/white cap Type A	250ml	PP
PL1521-2		500ml	PP
PL1521-3		1000ml	PP
PL1522-1	Washing bottle with red/white cap Type B	500ml	PP
PL1522-2		1000ml	PP
PL1523-1	Bottle with small mouth	250ml	PP
PL1523-2		500ml	PP
PL1523-3		1000ml	PP
PL1524-1	Bottle with big mouth, with gear cap	650ml	PP
PL1524-2		1000ml	PP
PL1524-3		1200ml	PP
PL1525-1	Bottle with big mouth, with round cap	30ml	PP
PL1525-2		50ml	PP
PL1525-3		100ml	PP
PL1525-4		200ml	PP
PL1525-5		250ml	PP

Cod.	Specification	Capacity	Material
PL1525-6	Bottle with big mouth, with round cap	300ml	PP
PL1525-7		500ml	PP
PL1525-8		620ml	PP
PL1525-9		1000ml	PP
PL1525-10		2000ml	PP
PL1526-1	Reagent Bottle	100ml	PP
PL1526-2		250ml	PP
PL1526-3		500ml	PP
PL1526-4	Media Bottle	1000ml	PP
PL1527-1		125ml	PET
PL1527-2		250ml	PET
PL1527-3	Funnel, short stem	500ml	PET
PL1527-4		1000ml	PET
PL1528-1		φ60x68mm	PP
PL1528-2	Funnel, long stem	φ67.5x74mm	PP
PL1529-1		φ60mm	PP
PL1529-2		φ75mm	PP
PL1529-3		φ90mm	PP

Reagent Bottle

- Made from medical PP and LDPE
- Various of volume and color for choices
- Different color of cap for choice: white, red, blue, green, yellow
- Excellent sealing design is leak-proof



Cod.	Description	Specification	Material	Body Color	Packing Info.
PL1601	4ml	Ø13 x 33mm	PE, PP	Natural / Black	5000pcs/ctn
PL1602	6ml	Ø14 x 40mm	PE, PP	Natural / Black	4000pcs/ctn
PL1603	10ml	Ø25 x 49mm	PE, PP	Natural / Black	3000pcs/ctn
PL1604	20ml	Ø27 x 62mm	PE, PP	Natural	2000pcs/ctn
PL1605	30ml	Ø35 x 63mm	PE, PP	Natural	1000pcs/ctn
PL1606	50ml	Ø43 x 65mm	PE, PP	Natural	1000pcs/ctn
PL1607	100ml	Ø43 x 114mm	PE, PP	Natural	500pcs/ctn
PL1608	15ml	Ø30.7 x 55mm	PP	Transparent	2000pcs/ctn
PL1609	100ml	Ø39.8 x 115mm	PP	Transparent	500pcs/ctn
PL1610	15ml	Ø24.5 x 62mm	PE, PP	Natural	1000pcs/ctn
PL1611	30ml	Ø33 x 65mm	PE, PP	Natural	1000pcs/ctn
PL1612	60ml	Ø37.5 x 86mm	PE, PP	Natural	1000pcs/ctn
PL1613	15ml	Ø24.5 x 62mm	PE, PP	Transparent	1000pcs/ctn
PL1614	30ml	Ø33 x 65mm	PE, PP	Transparent	1000pcs/ctn
PL1615	60ml	Ø37.5 x 86mm	PE, PP	Transparent	1000pcs/ctn
PL1616	15ml	Ø24.5 x 62mm	PE, PP	Brown	1000pcs/ctn
PL1617	30ml	Ø33 x 65mm	PE, PP	Brown	1000pcs/ctn
PL1618	60ml	Ø37.5 x 86mm	PE, PP	Brown	1000pcs/ctn



Cod.	Volum	Specification	Material	Body Color
Screw Cap				
PL1619	6ml	Ø23 x 62.8mm	PP + PE	Orange / Natural
PL1620	8ml	Ø26 x 66mm	PP + PE	Orange / Natural
PL1621	13ml	Ø38 x 65mm	PP + PE	Orange / Natural
PL1622	16ml	Ø38 x 80mm	PP + PE	Orange / Natural
PL1623	20ml	Ø45 x 57.5mm	PP + PE	Orange / Natural
PL1624	30ml	Ø45 x 65.5mm	PP + PE	Orange / Natural
PL1625	40ml	Ø45 x 96.5mm	PP + PE	Orange / Natural
PL1626	60ml	Ø45 x 140mm	PP + PE	Orange / Natural

Cod.	Volum	Specification	Material	Body Color
Pop Top Vial				
PL1627	20ml	Ø32.5 x 70mm	PP	Orange / Natural
PL1628	50ml	Ø39 x 70mm	PP	Orange / Natural
PL1629	80ml	Ø49 x 70mm	PP	Orange / Natural
PL1630	100ml	Ø49 x 135mm	PP	Orange / Natural




Cod.	Volum	Specification	Drop	Material	Body Color
PL1631	3ml	Ø14 x 47.5mm	40ul	PE, PP	Transparent
PL1632	3ml	Ø14 x 47.5mm	40ul	PE, PP	White
PL1633-1	6ml	Ø17.5 x 55mm	25ul	PE, PP	Transparent
PL1633-2	6ml	Ø17.5 x 55mm	40ul	PE, PP	Transparent
PL1633-3	6ml	Ø17.5 x 55mm	50ul	PE, PP	Transparent
PL1634-1	6ml	Ø17.5 x 55mm	25ul	PE, PP	White

Cod.	Volum	Specification	Drop	Material	Body Color
PL1634-2	6ml	Ø17.5 x 55mm	40ul	PE, PP	White
PL1634-3	6ml	Ø17.5 x 55mm	50ul	PE, PP	White
PL1635-1	7ml	Ø21 x 55mm	25ul	PE, PP	Transparent
PL1635-2	7ml	Ø21 x 55mm	40ul	PE, PP	Transparent
PL1635-3	7ml	Ø21 x 55mm	50ul	PE, PP	Transparent



Cod.	Volum	Specification	Drop	Material	Drop Vol.
PL1636	10ml	Ø24 x 69.4mm	40ul	PE, PP	Transparent
PL1637	10ml	Ø24 x 69.4mm	40ul	PE, PP	White
PL1636-1	10ml	Ø24 x 69.4mm	28ul	PE, PP	Transparent
PL1638	10ml	Ø25.2 x 56.6mm	70ul	PE, PP	White
PL1639	10ml	Ø25.2 x 56.6mm	70ul	PE, PP	Transparent
PL1640	20ml	Ø25 x 79mm	70ul	PE, PP	White
Dropper Bottle with Safety Cap					
PL1641	5ml	Ø12.5 x 55mm	50ul	PE, PP	Transparent
PL1642	8ml	Ø12.5 x 60mm	50ul	PE, PP	Transparent
PL1643	10ml	Ø12.5 x 66mm	50ul	PE, PP	Transparent
PL1644	5ml	Ø12.5 x 55mm	50ul	PE, PP	White
PL1645	8ml	Ø12.5 x 60mm	50ul	PE, PP	White
PL1646	10ml	Ø12.5 x 66mm	50ul	PE, PP	White

Cod.	Volum	Specification	Drop	Material	Drop Vol.
Push Type Dropper Bottle					
PL1647	5ml	Ø19 x 56mm	50ul	PE, PP	Transparent
PL1648	8ml	Ø20 x 61mm	50ul	PE, PP	Transparent
PL1649	10ml	Ø22 x 67mm	50ul	PE, PP	Transparent
PL1650	15ml	Ø25 x 70mm	50ul	PE, PP	Transparent
PL1651	5ml	Ø19 x 57mm	50ul	PE, PP	White
PL1652	8ml	Ø20 x 62mm	50ul	PE, PP	White
PL1653	10ml	Ø22 x 68mm	50ul	PE, PP	White
PL1654	15ml	Ø25 x 71mm	50ul	PE, PP	White



Cod.	Cap Color	Specification	Material	Body Color
PL1655-1	Green	Ø15.2 x 75.8mm	PE, PP, ABS	Transparent
PL1655-2	Blue	Ø15.2 x 75.8mm	PE, PP, ABS	Transparent
PL1655-3	Red	Ø15.2 x 75.8mm	PE, PP, ABS	Transparent
PL1655-4	Purple	Ø15.2 x 75.8mm	PE, PP, ABS	Transparent
PL1656	Green	Ø15.2 x 70mm	PE, PP, ABS	Transparent
PL1657	Purple	Ø15.2 x 75.8mm	PE, PP, ABS	Transparent


Cod.	Cap Color	Specification	Material	Body Color
PL1658-1	Blue	Ø17 x 110mm	PE, PP, ABS	Transparent
PL1658-2	Red	Ø17 x 110mm	PE, PP, ABS	Transparent
PL1659-1	Yellow	Ø12.4 x 66mm	PE, PP, ABS	Transparent
PL1659-2	Red	Ø12.4 x 66mm	PE, PP, ABS	Transparent
PL1660	White	Ø17 x 100mm	PE, PP, ABS	Transparent

Round Bottles, Narrow-Mouth



Cod.	Capacity	Closure Size	Material
PL1661	8ml	20mm	HDPE
PL1662	15ml	20mm	HDPE
PL1663	25ml	20mm	HDPE
PL1664	30ml	20mm	HDPE
PL1665	50ml	20mm	HDPE
PL1666	60ml	20mm	HDPE
PL1667	100ml	24mm	HDPE
PL1668	125ml	24mm	HDPE
PL1669	150ml	24mm	HDPE

Cod.	Capacity	Closure Size	Material
PL1670	200ml	24mm	HDPE
PL1671	250ml	24mm	HDPE
PL1672	500ml	28mm	HDPE
PL1673	1000ml	38mm	HDPE
PL1674	50ml	--	PTFE
PL1675	100ml	--	PTFE
PL1676	250ml	--	PTFE
PL1677	500ml	--	PTFE
PL1678	1000ml	--	PTFE



Cod.	Capacity	Closure Size	Material
PL1679	15ml	28mm	HDPE
PL1680	25ml	28mm	HDPE
PL1681	30ml	28mm	HDPE
PL1682	50ml	28mm	HDPE
PL1683	60ml	28mm	HDPE
PL1684	100ml	38mm	HDPE
PL1685	125ml	38mm	HDPE
PL1686	150ml	38mm	HDPE
PL1687	200ml	43mm	HDPE

Cod.	Capacity	Closure Size	Material
PL1688	250ml	43mm	HDPE
PL1689	500ml	53mm	HDPE
PL1690	1000ml	63mm	HDPE
PL1691	50ml	--	PTFE
PL1692	100ml	--	PTFE
PL1693	250ml	--	PTFE
PL1694	500ml	--	PTFE
PL1695	1000ml	--	PTFE



Cod.	Capacity	Closure Size	Material
Square Bottles, Narrow-Mouth			
PL1696	10ml	20mm	PE, PP
PL1697	25ml	20mm	PE, PP
PL1698	30ml	20mm	PE, PP
PL1699	50ml	24mm	PE, PP
PL16100	60ml	24mm	PE, PP
PL16101	100ml	38mm	PE, PP
PL16102	125ml	38mm	PE, PP
PL16103	250ml	38mm	PE, PP
PL16104	500ml	38mm	PE, PP
PL16105	1000ml	38mm	PE, PP

Cod.	Capacity	Closure Size	Material
Square Bottles, Wide-Mouth			
PL16106	60ml	28mm	PE, PP
PL16107	175ml	38mm	PE, PP
PL16108	250ml	43mm	PE, PP
PL16109	500ml	53mm	PE, PP
PL16110	1000ml	63mm	PE, PP
Carboys with Spigot			
PL16111	5L	38mm	PE, PP
PL16112	10L	38mm	PE, PP
PL16113	20L	83mm	PE, PP
PL16114	25L	83mm	PE, PP




Cod.	Capacity	Closure Size	Material
Rectangular Bottles			
PL16115	125ml	28mm	HDPE
PL16116	175ml	38mm	HDPE
PL16117	250ml	38mm	HDPE
PL16118	500ml	48mm	HDPE
PL16119	1000ml	53mm	HDPE
PL16120	2000ml	63mm	HDPE
Reagent Bottles with Stoppers, Narrow-Mouth			
PL16121	100ml	14/23mm	PP
PL16122	250ml	19/26mm	PP
PL16123	500ml	24/29mm	PP

Cod.	Capacity	Closure Size	Material
Reagent Bottles with Stoppers, Narrow-Mouth			
PL16124	1000ml	29/32mm	PP
PL16125	2000ml	29/32mm	PP
Reagent Bottles with Stoppers, Wide-Mouth			
PL16126	100ml	29/32mm	PP
PL16127	250ml	34/35mm	PP
PL16128	500ml	40/45mm	PP
PL16129	1000ml	60/64mm	PP
PL16130	2000ml	61/64mm	PP

Glassware

GG17 Beaker

- Low form, with spout, with printed graduations for 50-5000ml



Cod.	Capacity(ml)	Height (mm)	O.D. (mm)
PL1701	50ml	58mm	46mm
PL1702	100ml	72mm	52mm
PL1703	150ml	80mm	57mm
PL1704	200ml	88mm	64mm
PL1705	250ml	94mm	69mm
PL1706	300ml	102mm	76mm
PL1707	400ml	110mm	82mm
PL1708	500ml	115mm	87mm
PL1709	600ml	128mm	95mm
PL1710	800ml	138mm	102mm
PL1711	1000ml	150mm	110mm
PL1712	2000ml	200mm	113mm
PL1713	3000ml	210mm	150mm
PL1714	5000ml	280mm	178mm


GG17 Conical Flask

- Narrow neck



Cod.	Capacity	Height	O.D. of bulb	O.D. of neck
PL1726	50	85	51	22
PL1727	100	105	64	22
PL1728	150	120	71	28
PL1729	200	132	75	26
PL1730	250	145	82	30
PL1731	300	155	90	30
PL1732	500	195	100	35
PL1733	1000	225	128	40
PL1734	2000	285	170	53
PL1735	3000	330	195	54
PL1736	5000	380	235	65

Petri Culture Dish



Cod.	Bottom dia.(mm)	Cover dia.(mm)	Height (mm)
PL1747	60	65	15
PL1748	75	82	15
PL1749	90	98	18
PL1750	100	107	20
PL1751	120	130	25
PL1752	150	160	30

GG17 Boiling Flask


- Flat bottom, long narrow neck



Cod.	Capacity	Height	O.D. of bulb	O.D. of neck
PL1715	50ml	100mm	53	20
PL1716	100ml	120mm	65	21
PL1717	150ml	140mm	76	22
PL1718	200ml	150mm	83	23
PL1719	250ml	160mm	88	25
PL1720	300ml	170mm	92	26
PL1721	500ml	200mm	110	30
PL1722	1000ml	250mm	140	35
PL1723	2000ml	300mm	175	40
PL1724	3000ml	370mm	199	50
PL1725	5000ml	400mm	235	55

Condenser

- Allihn with bulbed inner tube



Cod.	Length of jacket(mm)	O.D. of upper tube(mm)	Cooling area(dm ²)	Length(mm)
PL1736	200	18	1	350
PL1738	250	18	1.5	420
PL1739	300	20	1.8	480
PL1740	400	22	2.7	590
PL1741	500	24	3.7	710
PL1742	600	24	4.9	800
PL1743	800	28	7.6	1050
PL1744	1000	28	10	1250
PL1745	1200	33	12.7	1360
PL1746	1500	33	17	1720

Vacuum Desiccator

- Clear, with ground-in stopcock and porcelain plate



Cod.	I.D.(mm)	Depth to plate support(mm)	Total height(mm)	Diameter of porcelain plate(mm)
PL1753	150	80	280	135
PL1754	180	100	350	150
PL1755	210	120	390	185
PL1756	240	130	440	210
PL1757	300	165	550	275

Glassware

Reagent Bottle


- Clear, narrow mouth, with ground-in glass stopper, or plastic stopper



Cod.	Capacity(ml)	Dia. Of bottle(mm)	Dia. Of mouth(mm)	Height(mm)
PL1758	30	40	18	76
PL1759	60	46	22	85
PL1760	125	57	24	110
PL1761	250	70	27	135
PL1762	500	85	33	172
PL1763	1000	106	38	202
PL1764	2500	145	48	270
PL1765	5000	185	58	338
PL1766	10000	225	68	430
PL1767	20000	265	82	520

Measuring Cylinder


- With spout



Cod.	Capacity (ml)	Graduation (ml)	Capacity tolerance(±ml)	Height (mm)
PL1768	5	0.1	0.1	110
PL1769	10	0.2	0.2	135
PL1770	25	0.5	0.5	160
PL1771	50	1	0.5	195
PL1772	100	1	1	250
PL1773	250	2or5	2	300
PL1774	500	5	5	350
PL1775	1000	10	10	430
PL1776	2000	20	20	500

Separatory Funnel


- Squib pear shape, with ground-in glass stopper and stopcock



Cod.	Capacity(ml)	O.D. of body(mm)	Total length(mm)	Diameter(mm)	I.D. of neck(mm)
PL1777	60	48	205	7	14
PL1778	125	67	260	7	19
PL1779	250	74	310	8	19
PL1780	500	90	355	9	24
PL1781	1000	120	430	10	29

Measuring Pipette

- With coding band




Cod.	Capacity (ml)	Graduation (ml)	Capacity tolerance(±ml)	Flow rate(s)	Height (mm)
PL1782	0.1	0.001	0.004	2-5	320
PL1783	0.2	0.002	0.006	2-5	320
PL1784	0.25	0.0025	0.008	2-5	320
PL1785	0.5	0.005	0.01	2-5	320
PL1786	1	0.01	0.015	3-6	350
PL1787	2	0.02	0.025	3-6	350
PL1788	5	0.05	0.05	5-10	350
PL1789	10	0.1	0.1	5-10	350
PL1790	25	0.1	0.1	11-21	470
PL1791	50	0.2	0.2	15-25	480

Reagent Bottle



Cod.	Capacity (ml)	Mouth Dia. (mm)	Spec.(mm)
PL1792	100	30	Ø56 x 100
PL1793	250	30	Ø70 x 138
PL1794	500	30	Ø86 x 176
PL1795	1000	30	Ø101 x 225




Cod.	Model	Volume Range	Accuracy	Repeatability
Number Reading Variable Volume Pipette				
PL1801	A2	0.1-2µl	±12.0-1.5%	≤6.00-0.70%
PL1802	A10	0.5-10µl	±2.5-1.0%	≤1.50-0.40%
PL1803	A20	2-20µl	±2.5-1.0%	≤1.50-0.30%
PL1804	A50	5-50µl	±2.0-1.0%	≤1.50-0.25%
PL1805	A100	10-100µl	±1.8-0.8%	≤0.50-0.15%
PL1806	A200	20-200µl	±1.8-0.6%	≤0.50-0.15%
PL1807	A1000	100-1000µl	±1.5-0.6%	≤0.30-0.15%
PL1808	A5000	500-5000µl	±1.2-0.6%	≤0.30-0.15%
PL1809	A10ml	1-10ml	±2.5-0.6%	≤0.30-0.15%
PL1810	AF5	5µl	±2.0%	≤0.80%
PL1811	AF10	10µl	±1.0%	≤0.50%
PL1812	AF20	20µl	±1.0%	≤0.30%
PL1813	AF25	25µl	±1.0%	≤0.30%
PL1814	AF50	50µl	±0.6%	≤0.30%
PL1815	AF100	100µl	±0.6%	≤0.25%
PL1816	AF200	200µl	±0.6%	≤0.15%
PL1817	AF500	500µl	±0.6%	≤0.20%
PL1818	AF1000	1000µl	±0.6%	≤0.15%
Direct Reading Variable Volume Pipette				
PL1819	H2	0.1-2µl	±12.0-1.5%	≤6.00-0.70%
PL1820	H10	0.5-10µl	±2.5-1.0%	≤1.50-0.40%
PL1821	H20	2-20µl	±2.5-1.0%	≤1.50-0.30%
PL1822	H50	5-50µl	±2.0-1.0%	≤1.50-0.25%
PL1823	H100	10-100µl	±1.8-0.8%	≤0.50-0.15%
PL1824	H200	20-200µl	±1.8-0.6%	≤0.50-0.15%
PL1825	H1000	100-1000µl	±1.5-0.6%	≤0.30-0.15%
PL1826	H5000	1000-5000µl	±1.2-0.6%	≤0.30-0.15%
PL1827	H10ML	1-10ml	±2.5-0.6%	≤0.30-0.15%
Fixed Volume Pipette				
PL1828	G5	5µl	±2.0%	≤0.80%
PL1829	G10	10µl	±1.0%	≤0.50%
PL1830	G20	20µl	±1.0%	≤0.30%
PL1831	G25	25µl	±1.0%	≤0.30%
PL1832	G50	50µl	±0.8%	≤0.30%
PL1833	G100	100µl	±0.8%	≤0.25%
PL1834	G200	200µl	±0.8%	≤0.15%
PL1835	G500	500µl	±0.8%	≤0.20%
PL1836	G1000	1000µl	±0.8%	≤0.15%

Cod.	Model	Volume Range	Accuracy	Repeatability
Number Reading Variable Volume Pipette				
PL1837	P2	0.1-2µl	±12.0-1.5%	≤6.00-0.70%
PL1838	P10	1-10µl	±2.5-1.0%	≤1.50-0.40%
PL1839	P20	2-20µl	±2.5-1.0%	≤1.50-0.30%
PL1840	P50	5-50µl	±2.0-1.0%	≤1.50-0.25%
PL1841	P100	10-100µl	±1.8-0.8%	≤0.50-0.15%
PL1842	P200	20-200µl	±1.8-0.6%	≤0.50-0.15%
PL1843	P1000	100-1000µl	±1.5-0.6%	≤0.30-0.15%
PL1844	P5000	1000-5000µl	±1.2-0.6%	≤0.30-0.16%
PL1845	P10ML	1-10ml	±2.5-0.6%	≤0.30-0.17%
PL1846	PF5	5µl	±2.0%	≤0.80%
PL1847	PF10	10µl	±1.0%	≤0.50%
PL1848	PF20	20µl	±1.0%	≤0.30%
PL1849	PF25	25µl	±1.0%	≤0.30%
PL1850	PF50	50µl	±0.6%	≤0.30%
PL1851	PF100	100µl	±0.6%	≤0.25%
PL1852	PF200	200µl	±0.6%	≤0.15%
PL1853	PF500	500µl	±0.6%	≤0.20%
PL1854	PF1000	1000µl	±0.6%	≤0.15%
Number Reading Variable Volume Pipette				
PL1855	U2	0.2-2µl	±12.0-2.0%	≤6.00-0.70%
PL1856	U10	1-10µl	±3.0-1.0%	≤1.50-0.40%
PL1857	U20	2-20µl	±3.0-1.0%	≤1.50-0.30%
PL1858	U50	5-50µl	±2.0-1.0%	≤1.50-0.20%
PL1859	U100	10-100µl	±2.0-0.8%	≤0.50-0.15%
PL1860	U200	20-200µl	±2.0-0.8%	≤0.50-0.15%
PL1861	U1000	100-1000µl	±1.5-0.8%	≤0.30-0.15%
PL1862	U5000	1000-5000µl	±1.0-0.7%	≤0.30-0.15%
PL1863	U10ML	1-10ml	±3.0-0.8%	≤0.30-0.15%
Majority Thoroughfare Variable Pipette				
PL1864	P8x20	2-20µl	±6.0-1.5%	≤2.50-0.40%
PL1865	P8x50	5-50µl	±3.0-1.2%	≤1.50-0.40%
PL1866	P8x200	20-200µl	±2.5-1.0%	≤1.25-0.50%
PL1867	P12x20	2-20µl	±6.0-1.5%	≤1.25-0.51%
PL1868	P12x50	5-50µl	±3.0-1.2%	≤1.25-0.52%
PL1869	P12x200	20-200µl	±2.5-1.0%	≤1.25-0.53%
PL1870	A8x20	2-20µl	±6.0-1.5%	≤2.50-0.40%
PL1871	A8x50	5-50µl	±3.0-1.2%	≤1.50-0.40%
PL1872	A8x200	20-200µl	±2.5-1.0%	≤1.25-0.50%
PL1873	A12x20	2-20µl	±6.0-1.5%	≤1.25-0.51%
PL1874	A12x50	5-50µl	±3.0-1.2%	≤1.25-0.52%
PL1875	A12x200	20-200µl	±2.5-1.0%	≤1.25-0.53%
Repetitive Pipette				
HX Tips				
PL1876	0.5ml	10-50µl	±0.9%	≤1.4% to ≤0.6%
PL1877	1.25ml	25-125µl	±0.9%	≤0.8% to ≤0.3%
PL1878	2.5ml	50-250µl	±0.8%	≤0.8% to ≤0.2%
PL1879	5.0ml	100-500µl	±0.6%	≤0.6% to ≤0.2%
PL1880	12.5ml	250-1250µl	±0.4%	≤0.4% to ≤0.2%
PL1881	25ml	500-2500µl	±0.4%	≤0.4% to ≤0.2%
PL1882	50ml	1000-5000µl	±0.4%	≤0.4% to ≤0.15%

Electromotion Pipette



Cod.	Description
PL1883	Suitable for 0.1ml-100ml large range liquid operation European style charger 100-240V/50HZ Once charge up, continuous working for 8 hours



Cod.	Description	Color
Pipette Pump		
PL1884	2ml	purple
PL1885	10ml	green
PL1886	25ml	red
Large Volume Pipette Controller		
PL1887	0.1ml-100ml	grey
PL1888	0.1ml-100ml	blue
PL1889	0.1ml-100ml	red
PL1890	0.1ml-100ml	green
PL1891	0.1ml-100ml	yellow
Finland Pipette		
PL1892	0.5-10µl	/
PL1893	5-50µl	/
PL1894	20-200µl	/
PL1895	100-1000µl	/
PL1896	1000-5000µl	/




Cod.	Description
PL1897	Full/Half autoclavable mechanical micropipette pipettes



Cod.	Description
PL1898	Full/Half autoclavable mechanical micropipette multi-channel pipettes

Pipette Rack



Cod.	Description	Specification
PL1899	Pipette rack	9 pipettes
PL18100		9 pipettes
PL18101		5 pipettes

Cod.	Description	Specification
PL18102	Pipette rack	5 pipettes
PL18103	Multi-use rack	Suitable for pipette, serological pipette and ESR pipette

Instruments

JSQA Hemocytometer

The Instrument is microcomputer controlled hemocytometer, which could be counted, numerical statement and display. It is 10 groups of 2 digital. It could display ten groups of numbers and totals accurately, satisfy the clinic to test the blood routine detection. It is attractive appearance, rational construction, clear data, it is convenient to use, especially applicable to the middle size and small size hospitals.




J-3 Colonometer

This Bacteria Colony Counter is made of CMOS integrated electric, is a digital display for counting the germs, capacity of the counter is from 0 to 999, the button return back to "0" automatically. When error in counter, push the "-" and "+", and the counter will be less or add. With sensor, it's sensitive, incl sound indicator, to ensure the counting is correct. This instrument supply with special pen, enhance the efficiency of the counter, the magnifier can be adjusted to perfect position with flexible arm. The magnification up to 3x, it's easy to compare bacteria, apply to Petri dish from 50Ø~150Ømm.



WJX-1 Hemoglobinometer

The product is assembled with the micro-computer with the characters of intelligence, inletting adjustable, easy for operation and accurate testing as well as the function stability. It is an ideal successor to the simple old hemoglobinometers.




PSMC-I Decapper for Vacuum Blood Collection Tubes

It would be simpler and more convenient by using the automatic machine to reduce medical staff's labor intensity and improve the work efficiency. Effectively avoid being infected by skin scratch when the staff open the cap of cracked glass tube by hand. Due to the negative pressure inside the tube, too much force during the cap opening will cause blood spatter. The machine successfully solves the danger of cross infection and avoid the trouble of frequent hand washing. Radically preclude the danger of potential biological contamination caused by residual. Radically preclude the danger of danger of potential biological contamination caused by residual inside the tube during manual operation. Meanwhile, it is highly meaningful to bio-safety.



PSMR-III Multifunction Mixer

Multi-Mixer has the capability to perform three modes of mixing, all of which are programmable. It is very simple to operate within seconds. It can mix variable substance such as blood sample, urine sample, injectant, blood plasma etc.



PSMR-II Blood Roller Mixer

The instrument is used for complete mixing of blood, powder reagents, etc. in the anticoagulation tubes or vials. It is proved to be more accurate and convenient than manual mixing by Shanghai Center for Clinical Laboratory as well as the big hospital all over China. According to different specimen, we can provide various series, such as rolling and swing, CW/CCW rotation and speed adjustable machine. Also, the hollow shaft can be lengthened with limits to increase the mixing capacity. It is an ideal instrument for clinical lab.



Instruments


PSMR-IV Blood Mixer

Equipped with silicagel, it's easy and convenient for putting out specimen tubes. The working plate runs circularly, speed can be adjusted by the operating staff according to different specimen. It is applicable to the mixing of clinical blood routine and vacuum blood collection tubes.



PSMR-V 3D Rotation Mixer

The kind of Mixer with three-dimensional motion provides gentle and thorough mixing. It is ideal for sample mixing in centrifuge tubes, as well as re-suspension of large, fragile DNA FRAGMENTS. The HXMR-V accommodates most size tube, centrifuge tube, microwell plate, petri dishes and tissue culture plate and cell culture bottle. This equipment is suitable to work in the incubator or cold room.



XH-C Vortex

Being of the compact design, no noise and easy for operation, the product has two functions of continuous and inching, which is suitable for blending of the blood routine (anticoagulation tube, 820 sample cup and 900 sample cup etc.) with character of full mixing and no damage to the blood cells. Comparing to the manual operation, it is more convenient and also can be used for liquid mixing in other containers.



PS-201BY Rocking Mixer

It's a machine which swings up and down with little swing amplitude. Frequency adjustable. It's used for fixing of running gel spectral band, rocking of Coomassie blue and decolorization and fixing of silver nitrate staining, etc. It's used for development and fixation of X-ray film, further treatment cellulose membrane after electrophoretic migration, molecular hybridization, cell culture and blood mixing.



XH-B Vortex

The product is equipped with imported motor, which is used for mixing of the liquids in centrifuge tubes, test tubes and breakers. According to the operator's strength, the vortex frequency is adjustable to ensure the full and fast mixing. The working plate is designed to be hard wearing and easy for cleaning.



XH-D Vortex

The economical vortex mixer is in a wide range of operation, there are two kinds of working plate: rubber and plastic. It can be used for liquids mixing in different test tubes, flasks, beakers. Dual mode switch makes the vortex work continuously or in inching. The speed of both models could be adjusted according to different requirement. High efficiency and low noise.



Instruments

PS-201A Oscillator (Minim-Oscillator)

Use the silicone tray, which can be used for different culture plates, micro-plate oscillating simultaneously with high effectiveness and speed adjustable. It is the successor to the old products.



PS-201C Oscillator

Professional design and no noise, 6 plates can be placed in at one times, speed adjustable. Imported motor, little noise, could equipped with digital display of speed and timing if needed.



PS-201BD Oscillator (horizontal rotor)

The product is used for mixing of RPR testing card with character of stability, no noise, reasonable rotating speed and compact design. It is applicable to the epidemic prevention station, blood station, dermatosis clinic and hospital laboratories.



PS-202 Oscillator (Powder mixer)

It is professional design, there is no shift and noise when it is working, 19 different kinds of bottles can be put in at one time and easy for operation. It can mix the powder medicine clearly within 20 seconds which is hard to dissolve by the normal way, avoid the waste of the medicine. It's a good helper for the staff in hospitals; it could be used for mixing of samples of liquid sample also, the number and diameter of holes could be customized, to adapt different container's mix. Time fixed, speed adjustable, could work 12 hours continuously, it is convenient to improve work efficiency.



PS-201BS Oscillator (horizontal rotor)

The product is used for mixing of RPR testing card with character of stability, no noise, reasonable rotating speed and compact design. It is applicable to the epidemic prevention station, blood station, dermatosis clinic and hospital laboratories.



DSZ-II Digital Timing Clock

15 groups timing at the same time, 12 pieces of melody for choice; any one second could be forecasted; display the data and time; with features of memory protection after power failure; high definition, high brightness VALCD display.



Vacuum tube

- Safety cap
- High quality tube
- High efficiency
- High quality additives
- Customize the label



No Additive Tube



Pro-coagulation Tube



EDTA Tube



Heparin Tube



ESR Tube



Glucose Tube



Gel and Clot Activator Tube




PT Tube

Cod.	Specification	Volume	Additive	Color	Tube Material
	No Additive Tube				
PM0101	Ø13x75mm	3ml	No	Red	Glass / PET
PM0102	Ø13x75mm	4ml			Glass / PET
PM0103	Ø13x75mm	5ml			Glass / PET
PM0104	Ø13x100mm	5ml			Glass / PET
PM0105	Ø13x100mm	6ml			Glass / PET
PM0106	Ø13x100mm	7ml			Glass / PET
PM0107	Ø16x100mm	8ml			Glass / PET
PM0108	Ø16x100mm	9ml			Glass / PET
PM0109	Ø16x100mm	10ml			Glass / PET
	Pro-coagulation Tube				
PM0110	Ø13x75mm	3ml	Clot Activator	Orange	Glass / PET
PM0111	Ø13x75mm	4ml			Glass / PET
PM0112	Ø13x75mm	5ml			Glass / PET
PM0113	Ø13x100mm	5ml			Glass / PET
PM0114	Ø13x100mm	6ml			Glass / PET
PM0115	Ø13x100mm	7ml			Glass / PET
PM0116	Ø16x100mm	8ml			Glass / PET
PM0117	Ø16x100mm	9ml			Glass / PET
PM0118	Ø16x100mm	10ml			Glass / PET
	EDTA Tube				
PM0119	Ø13x75mm	1ml	EDTA K2	Light Purple	Glass / PET
PM0120	Ø13x75mm	2ml			Glass / PET
PM0121	Ø13x75mm	3ml			Glass / PET
PM0122	Ø13x75mm	4ml			Glass / PET
PM0123	Ø13x75mm	5ml			Glass / PET
PM0124	Ø13x100mm	5ml			Glass / PET
PM0125	Ø13x100mm	6ml			Glass / PET
PM0126	Ø13x100mm	7ml			Glass / PET
PM0127	Ø16x100mm	10ml			Glass / PET
PM0128	Ø13x75mm	1ml	EDTA K3	Purple	Glass / PET
PM0129	Ø13x75mm	2ml			Glass / PET
PM0130	Ø13x75mm	3ml			Glass / PET
PM0131	Ø13x75mm	4ml			Glass / PET
PM0132	Ø13x75mm	5ml			Glass / PET
PM0133	Ø13x100mm	5ml			Glass / PET
PM0134	Ø13x100mm	6ml			Glass / PET
PM0135	Ø13x100mm	7ml			Glass / PET
PM0136	Ø16x100mm	10ml			Glass / PET
PM0137	Ø13x75mm	2ml	EDTA Na2	Deep Purple	Glass / PET
PM0138	Ø13x100mm	5ml			Glass / PET

Cod.	Specification	Volume	Additive	Color	Tube Material
	Heparin Tube				
PM0139	Ø13x75mm	3ml	Sodium Heparin	Green	Glass / PET
PM0140	Ø13x75mm	4ml			Glass / PET
PM0141	Ø13x75mm	5ml			Glass / PET
PM0142	Ø13x100mm	7ml			Glass / PET
PM0143	Ø16x100mm	10ml			Glass / PET
PM0144	Ø13x75mm	3ml	Lithium Heparin	Deep Green	Glass / PET
PM0145	Ø13x75mm	4ml			Glass / PET
PM0146	Ø13x75mm	5ml			Glass / PET
PM0147	Ø13x100mm	7ml			Glass / PET
PM0148	Ø16x100mm	10ml			Glass / PET
	ESR Tube				
PM0149	Ø13x75mm	1.6ml	3.8% Sodium Citrate	Black	Glass / PET
PM0150	Ø13x75mm	2.4ml			Glass / PET
PM0151	Ø9x120mm	1.28ml			Glass / PET
	Glucose Tube				
PM0152	Ø13x75mm	2ml	Sodium Fluoride + Sodium Heparin	Grey	Glass / PET
PM0153	Ø13x75mm	3ml			Glass / PET
PM0154	Ø13x75mm	4ml			Glass / PET
PM0155	Ø13x75mm	5ml			Glass / PET
PM0156	Ø13x100mm	3ml			Glass / PET
PM0157	Ø13x100mm	4ml			Glass / PET
PM0158	Ø13x100mm	5ml			Glass / PET
PM0159	Ø13x100mm	6ml			Glass / PET
	Gel and Clot Activator Tube				
PM0160	Ø13x75mm	3ml	Gel&Clot Activator	Yellow	Glass / PET
PM0161	Ø13x75mm	4ml			Glass / PET
PM0162	Ø13x100mm	5ml			Glass / PET
PM0163	Ø13x100mm	6ml			Glass / PET
PM0164	Ø16x100mm	7ml			Glass / PET
PM0165	Ø16x101mm	8ml			Glass / PET
PM0166	Ø16x102mm	9ml			Glass / PET
	PT Tube				
PM0167	Ø13x75mm	1.8ml	3.2% Sodium Citrate	Blue	Glass / PET
PM0168	Ø13x75mm	2.7ml			Glass / PET
PM0169	Ø13x75mm	3.6ml			Glass / PET
PM0170	Ø13x75mm	4.5ml			Glass / PET
PM0171	Ø13x75mm (Double deck)	1.8ml			Glass / PET
PM0172	Ø13x75mm (Double deck)	2.7ml		Glass / PET	

Non-Vacuum tube

- Made of medical grade PP material
- Siliconization treatment on the inner surface, no blood adhesion after collection
- High quality gel makes the serum separation complete with no fibrin adhering to the inner surface of the tube



Cod.	Specification	Additive	Color
PM0201	Ø13x75mm	EDTA.K2 / EDTA.K3	Purple
PM0202	Ø13x100mm		
PM0203	Ø13x75mm	Sodium Heparin / Lithium Heparin	Green
PM0204	Ø13x100mm		
PM0205	Ø13x75mm	Gel & Clot Activator	Yellow
PM0206	Ø13x100mm		

Cod.	Specification	Additive	Color
PM0207	Ø13x75mm	Sodium Fluoride+Heparin Sodium / Sodium Fluoride+EDTA.K2	Grey
PM0208	Ø13x100mm		
PM0209	Ø13x75mm	No	Red
PM0210	Ø13x100mm		

Micro Blood Collection Tube




Cod.	Vol.	Additive	Color
PM0211	0.5ml	EDTA.K2 / EDTA.K3	Pink
PM0212	0.5ml	Lithium Heparin / Sodium Heparin	Green
PM0213	0.5ml	No	Red
PM0214	0.5ml	Gel & Clot Activator	Yellow
PM0215	0.5ml	EDTA.K2	Purple
PM0216	0.5ml	Heparin Lithium / Heparin Sodium	Green
PM0217	0.5ml	Plain	Red
PM0218	0.5ml	Gel & Clot Activator	Yellow



Cod.	Description	Specification	Packing Info.
PM0219	Glass Capillary	blue coding, without heparin	500pcs/bottle
PM0220			100pcs/bottle
PM0221		red coding, with heparin	500pcs/bottle
PM0222			100pcs/bottle


Swab

- Round tip, compact and shred resistant, free of loose during sampling
- Supplied with sterile and non-sterile



Cod.	Description	Specification	Stick Material
PM0301	Transport Swab	Male	stainless steel
PM0302		Female	wooden
PM0303	Transport Swab	Tube without Media	plastic
PM0304			wooden
PM0305			stainless steel
PM0306			bamboo

Cod.	Description	Specification	Stick Material
PM0307	Collection Swab	Ø2 x 25mm	ABS + Nylon Flocked
PM0308		Ø2.5 x 18mm	
PM0309		Ø5 x 18mm	
PM0310		Ø5 x 25mm	
PM0311		Ø6 x 28mm	PP + Nylon Flocked



Cod.	Description	Specification	Medium
PM0312	Transport Swab	Individual pack, Gamma sterilized	Amies Clr
PM0313			Amies
PM0314	Transport Swab	Individual paper-plastic pack, EO sterilized	Without Medium
PM0315			Amies
PM0315-1	Transport Swab	Individual paper-plastic pack, EO sterilized	Amies with Charcoal
PM0315-2			Stuart
PM0315-3			Gary Blair

Cod.	Description	Specification	Stick Material
PM0316	Cotton Swabs Absorbent Cotton	4"xø5mm	plastic
PM0317		4"xø6mm	wooden
PM0318	Swab packed with diayzing paper	individual pack, sterile	plastic or wooden

Disposable Cervical Scraper

- Ayres scraper, designed to collect cytology specimens from the cervix and is also adaptable to collect specimens from the tongue, throat or vulva
- Duo-End Scraper, enables the clinician to use either the Ayres Tips or modified Elongared Tip, providing cell samples from both cervix and endocervix in one step
- Elongated scraper, allows deeper penetration into the endocervical area




Cod.	Description	Material
PM0319	Type A	Plastic
PM0320	Type B	Plastic
PM0321	Type C	Plastic
PM0322	Type D	Plastic
PM0323	Type E	Plastic

Cod.	Description	Material
PM0324	Type F	Birch wood
PM0325	Type G	Birch wood
PM0326	Type H	Birch wood
PM0327	Type I	Birch wood
PM0328	Type J	Birch wood


Cervical Brush

- Various lengths of plastic hairs design, ideal for collecting cells from the entire required sampling area of the cervix, ectocervical and endocervical cells for Pap smears
- Plastic handle, with soft, flexible brush, providing more efficient sampling and reliable cell transfer



Cod.	Model	Description
PM0401	Type A-1	Push tube, pink brush, individually packaged, EO sterile.
PM0402	Type A-2	No push tube, pink brush, individually packaged, EO sterile.
PM0403	Type B-1	Push tube, blue brush, individually packaged, EO sterile.
PM0404	Type B-2	No push tube, blue brush, individually packaged, EO sterile.
PM0405	Type C	Blue brush bar, individually packaged, EO sterile.
PM0406	Type D	Blue brush bar, individually packaged, EO sterile.
PM0407	Type E	Fluffy head, blue brush bar, individually packaged, EO sterile.
PM0408	Type F	Fluffy head, blue brush bar, individually packaged, EO sterile.
PM0409	Type G	Blue brush bar, individually packaged, EO sterile.
PM0410	Type H	Fluffy head, blue brush bar, individually packaged, EO sterile.

Cod.	Model	Description
PM0411	Type I	Nylon head, blue brush bar, individually packaged, EO sterile.
PM0412	Type J	Nylon head, white brush bar, individually packaged, EO sterile.
PM0413	Type K	Two joints, white, flocking, individually packaged, EO sterile.
PM0414	Type L	Yellow brush bar, individually packaged, EO sterile.
PM0415	Type M	Nylon head with head white brush bar, individually packaged, EO sterile.
PM0416	Type N	Fluffy head, blue brush bar, individually packaged, EO sterile.
PM0417	Type O	White brush bar, individually packaged, EO sterile



Cod.	Model	Description
PM0418	Type P-1	Nylon head, blue brush bar, individually packaged, EO sterile
PM0419	Type P-2	Nylon head, blue brush bar, individually packaged, EO sterile
PM0420	Type P-3	Nylon head, white brush bar, individually packaged, EO sterile
PM0421	Type P-4	Nylon head with head white brush bar, individually packaged, EO sterile.
PM0422	Type P-5	Nylon head, white brush bar, individually packaged, EO sterile
PM0423	Type P-6	Nylon head with head white brush bar, individually packaged, EO sterile.
PM0424	Type P-7	Nylon head with head blue brush bar, individually packaged, EO sterile.
PM0425	Type P-8	Nylon head with head pink brush bar, individually packaged, EO sterile.
PM0426	Type P-9	Nylon head with head white brush bar, individually packaged, EO sterile.
PM0427	Type P-10	Nylon head with head white brush bar, individually packaged, EO sterile.
PM0428	Type Q-1	Nylon head, blue brush bar, individually packaged, EO sterile
PM0429	Type Q-2	Nylon head, white brush bar, individually packaged, EO sterile
PM0430	Type Q-3	Nylon head, green brush bar, individually packaged, EO sterile


Blood Lancet



Cod.	Description	Specification
PM0501	Safety Type Blood Lancet	21G
PM0502		23G
PM0503		26G
PM0504	Pressure Activated Blood Lancet	21G
PM0505		23G
PM0506		26G
PM0507	Twist Type Blood Lancet	---
PM0508	Pull Type Blood Lancet	---
PM0509	Flake Type Blood Lancet	---


Cod.	Description	Specification
PM0510	Blood Collection Needle	20G
PM0511		21G
PM0512		22G
PM0513	Soft link blood collection needle	23G
PM0514		21G
PM0515		22G
PM0516	Lancing Device	40ml
PM0517	Blood Collection Needle	18-23G
PM0518	Single-winged / Double-winged	
PM0519	Tourniquet with Buckle	--

Disposable Scalp Vein Sets



Cod.	Color	Size
PM0519	Grey	27G
PM0520	Brown	26G
PM0521	Orange	25G
PM0522	Purple	24G
PM0523	Blue	23G
PM0524	Black	22G
PM0525	Green	21G
PM0526	Yellow	20G
PM0527	Cream	19G
PM0528	Pink	18G

Disposable Hypodermic Needle



Cod.	Color	Size
PM0529	Grey	27G
PM0530	Brown	26G
PM0531	Orange	25G
PM0532	Purple	24G
PM0533	Blue	23G
PM0534	Black	22G
PM0535	Green	21G
PM0536	Yellow	20G
PM0537	Cream	19G
PM0538	Pink	18G

Syringe 3 Parts



Cod.	Description	Size
PM0601	Syringe Luer Slip, w/wo needle, 3 parts	1ml
PM0602		2ml
PM0603		3ml
PM0604		5ml
PM0605		10ml
PM0606		20ml
PM0607		30ml
PM0608		50ml
PM0609		60ml
PM0610		100ml

Syringe 3 Parts



Cod.	Description	Size
PM0611	Syringe Luer Lock, w/wo needle, 3 parts	1ml
PM0612		2ml
PM0613		3ml
PM0614		5ml
PM0615		10ml
PM0616		20ml
PM0617		30ml
PM0618		50ml
PM0619		60ml
PM0620		100ml

Syringe 2 Parts



Cod.	Description	Size
PM0621	Syringe Luer Slip, w/wo needle, 2 parts	1ml
PM0622		2ml
PM0623		3ml
PM0624		5ml
PM0625		10ml
PM0626		20ml
PM0627		30ml
PM0628		50ml
PM0629		60ml
PM0630		100ml

Syringe 2 Parts



Cod.	Description	Size
PM0631	Syringe Luer Lock, w/wo needle, 2 parts	1ml
PM0632		2ml
PM0633		3ml
PM0634		5ml
PM0635		10ml
PM0636		20ml
PM0637		30ml
PM0638		50ml
PM0639		60ml
PM0640		100ml

Safety AD Syringe



Cod.	Size	Spec.
PM0641	1ml	Luer Slip
PM0642	2ml	Luer Slip
PM0643		Luer lock
PM0644	3ml	Luer Slip
PM0645		Luer lock
PM0646	5ml	Luer Slip
PM0647		Luer lock
PM0648	10ml	Luer Slip
PM0649		Luer lock
PM0650	20ml	Luer Slip
PM0651		Luer lock

AD Syringe




Cod.	Size	Spec.
PM0652	0.5ml	Luer Slip
PM0653	1ml	Luer Slip
PM0654	2ml	Luer Slip
PM0655		Luer lock
PM0656	3ml	Luer Slip
PM0657		Luer lock
PM0658	5ml	Luer Slip
PM0659		Luer lock
PM0660	10ml	Luer Slip
PM0661		Luer lock
PM0662	20ml	Luer Slip
PM0663		Luer lock



Cod.	Description	Spec.
PM0664	Insulin Syringe Fixed Needle	0.3ml / 100ul
PM0665		0.5ml / 100ul
PM0666		1ml / 100ul

Cod.	Description	Spec.
PM0667	Insulin Syringe Fixed Needle	0.5ml / 40ul
PM0668		1ml / 40ul
PM0669	Ultra-safe Syringe for Animal	----



Cod.	Description	Spec.
PM0670	Disposable Infusion Sets	Type A
PM0671		Type B
PM0672		Type C
PM0673		Type D


Cod.	Description	Spec.
PM0674	Disposable Infusion Sets	Type E
PM0675		Type F
PM0676		Type G
PM0677		Type H




Cod.	Description	Spec.
PM0678	Disposable Infusion Sets	Type I
PM0679		Type J
PM0680		Type K
PM0681		Type L

Cod.	Description	Spec.
PM0682	Disposable Blood Transfusion Sets	Type A
PM0683		Type B
PM0684		Type C
PM0685		Type D

I.V. Catheter



Cod.	Description	Spec.	Size.
PM0686	Volumetric I-V Burette	Soft Chamber	100ml
PM0687		Soft Chamber	150ml
PM0688	Volumetric I-V Burette	Soft Chamber	100ml
PM0689		Hard Chamber	100ml




Cod.	Description	Color	Size
PM0690	Pen Type	Yellow	24G
PM0691		Blue	22G
PM0692		Pink	20G
PM0693		Green	18G
PM0694		Grey	16G
PM0695	Butterfly Type	Orange	14G
PM0696		Yellow	24G
PM0697		Blue	22G
PM0698		Pink	20G
PM0699		Green	18G
PM06100	Normal / Safety Type	Grey	16G
PM06101		Orange	14G
PM06102		Yellow	24G
PM06103		Blue	22G
PM06104		Pink	20G
PM06105		Green	18G
PM06106	Heparin Cap	--	--

Oxygen Mask

Swivel Nebulizer Mask

- Made of medical grade PVC
- 360° rotation connector
- Available with 7' anti-crush tubing.



Cod.	Size
PM0701	Adult
PM0702	Pediatric

Swivel Oxygen Mask


- Made of medical grade PVC
- 360° rotation connector
- Available with 7' anti-crush tubing.



Cod.	Size
PM0703	Adult
PM0704	Pediatric

Simple Oxygen Mask

- Made of medical PVC
- Adjustable nose clip
- Medium concentration
- Available with anti-crush 7ft tubing



Cod.	Size
PM0705	XL
PM0706	L
PM0707	M
PM0708	S

Nebulizer Mask

- Made of medical PVC
- Adjustable nose clip
- Available with anti-crush 7ft tubing



Cod.	Size
PM0709	XL
PM0710	L
PM0711	M
PM0712	S

Simple Oxygen Mask

- Made of medical PVC
- Without nose clip
- Available with anti-crush 7ft tubing
- Flexible connector, 360° rotation, more convenient for use



Cod.	Size
PM0713	Adult
PM0714	Pediatric

Nebulizer Kit

- Made of medical PVC, 100% latex free
- Available with four types nebulizers
- Available with anti-crush 7ft tubing
- Two colors for choice: Green and transparent



Cod.	Size
PM0715	T-type
PM0716	Y-type

Non Rebreathing Mask

- Made of medical PVC
- Adjustable nose clip
- Available with anti-crush 7ft tubing



Cod.	Size
PM0717	XL
PM0718	L
PM0719	M
PM0720	S

Adjustable Venturi Mask

- Made of medical PVC
- Adjustable nose clip
- Oxygen concentration is 24%-50%
- Available with anti-crush 7ft tubing



Cod.	Size
PM0721	XL
PM0722	L
PM0723	M
PM0724	S

Adjustable Venturi Mask

- Made of medical PVC
- Adjustable nose clip
- Venturi valve: blue 24%, yellow 28%, white 31%, green 35%, pink 40%, orange 50%
- Available with anti-crush 7ft tubing




Cod.	Size
PM0725	XL
PM0726	L
PM0727	M
PM0728	S

MEDICAL

Anesthesia Mask


- Made of medical PVC
- Adjustable nose clip



Cod.	Size
PM0729	Preemie
PM0730	Neonate
PM0731	Infant
PM0732	Toddler
PM0733	Child/Small adult
PM0734	Medium adult
PM0735	Large adult

Aerosol Mask


- Made of medical PVC
- Adjustable nose clip
- Available with nebulizer and oxygen supply tubing



Cod.	Size
PM0736	XL
PM0737	L
PM0738	M
PM0739	S

CO₂ Mask


- Simultaneous CO₂ monitoring and oxygen delivery
- Available with anti-crush 7ft tubing



Cod.	Size
PM0740	Adult
PM0741	Pediatric

CO₂ Cannula

- Tubing in peel design
- Tubing length: 7ft



Cod.	Size
PM0742	Adult
PM0743	Pediatric
PM0744	Infant

Soft-tip Nasal Oxygen Cannula

- Made of medical PVC
- Generally 7ft anti-crush PVC tubing
- Super soft tips to comfort the patient



Cod.	Size
PM0745	Adult
PM0746	Pediatric
PM0747	Infant
PM0748	Neonate

Double-Hole Nasal Cannula

- Made of medical PVC
- Two soft tips, over ear style
- Generally 200cm anti-crush PVC tubing



Cod.	Size
PM0749	adult
PM0750	pediatric
PM0751	small pediatric

Oxygen Connection Tubing


- Made of medical PVC
- Generally 7ft anti-crush PVC tubing



Cod.	Color
PM0752	transparent
PM0753	green

Disposable PVC Laryngeal Mask

- Made of medical PVC
- Single use only




Cod.	Size	Patient weight	Cuff vol.
PM0754	1.0	<5kgs	<4ml
PM0755	1.5	5-10kgs	<7ml
PM0756	2.0	10-20kgs	<10ml
PM0757	2.5	20-30kgs	<14ml
PM0758	3.0	30-50kgs	<20ml
PM0759	4.0	50-70kgs	<30ml
PM0760	5.0	70-100kgs	<40ml

Disposable Silicone Laryngeal Mask



Cod.	Size	Patient weight	Cuff vol.
PM0761	1.0	<5kgs	<4ml
PM0762	1.5	5-10kgs	<7ml
PM0763	2.0	10-20kgs	<10ml
PM0764	2.5	20-30kgs	<14ml
PM0765	3.0	30-50kgs	<20ml
PM0766	4.0	50-70kgs	<30ml
PM0767	5.0	70-100kgs	<40ml



Cod.	Size	LCD	Spec.
PM0801	127x18x10mm	15.5x6.5mm	/
PM0802	122x17x10mm	20x7.5mm	Large LCD display
PM0803	127x18x10mm	17.5x7.5mm	Transparent(blue, yellow,red)
PM0804	127x18x10mm	16.5x7mm	/
PM0805	133x18x13mm	20x7mm	Water-Proof, Large LCD Display
PM0806	125x17x10mm	16x6.5mm	Water-Proof
PM0807	127x18x10mm	17.5x7.5mm	Water-Proof
PM0808	127x18x10mm	15.5x6.5mm	High Accuracy

Cod.	Size	LCD	Spec.
PM0809	133x18x13mm	20x7mm	High Accuracy, Large LCD Display, Water-proof
PM0810	126x18x10mm	20x7mm	Flexible Probe, Large LCD Display
PM0811	126x18x10mm	20x7mm	Flexible Probe, Large LCD Display, Water-proof
PM0812	125x20x11mm	20x7mm	Flexible Probe, Large LCD Display, Water-proof



Cod.	Size	LCD	Spec.
PM0813	134x22x14mm	20.7x7mm	Fast Reading, Large LCD Display
PM0814	134x22x14mm	20.7x7mm	Fast Reading, Large LCD Display,Water-Proof
PM0815	134x22x14mm	20x7.5mm	Fast Reading, Large LCD Display,Water-Proof
PM0816	134x22x14mm	20.7x7mm	Fast Reading, Large LCD Display,Water-Proof
PM0817	134x22x14mm	20.7x7mm	Fast Reading, Large LCD Display,Flexible Probe,Water-Proof

Cod.	Size	LCD	Spec.
PM0818	134x22x14mm	20x7.5mm	Fast Reading, Large LCD Display,Flexible Probe,Water-Proof
PM0819	143x24x15.5mm	31x14mm	Water-proof, Flexible tip
PM0820	143x24x15.5mm	31x14mm	Water-proof, Flexible tip
PM0821	143x24x15.5mm	31x14mm	Water-proof, Flexible tip
PM0822	62x40x45mm	15x6mm	Baby Pacifier Thermometer
PM0823	62x46x50mm	15x6mm	Baby Pacifier Thermometer
PM0824	150.9x26x19mm	20x7mm	Character Digital Thermometer



Cod.	Description
	Infrared Ear Thermometer
PM0825	Reading in 1 second,10 memories,with 20 covers,need a new disposable probe cover for each use
PM0826	Reading in 1 second,10 memories,cover free design,easy-clean
PM0827	Reading in 1 second,10 memories,cover free design,easy-clean
PM0828	Reading in 1 second,10 memories,with dust cover,easy-clean

Cod.	Description
	Infra-red Forehead Thermometer
PM0829	Reading in 1 second, Ear / Forehead, 10 memories, need a new disposable probe cover for each use
PM0830	Reading in 1 second, 10 memories, quick and accurate
PM0831	Reading in 1 second, 32 memories, quick and accurate, Non-contact, Voice identification
PM0832	Reading in 1 second, 32 memories, quick and accurate, Non-contact, Voice identification, Self-set Fever-alarm

Vaginal Speculum

PM0901

Vaginal speculum with Smoke extractor type
(For gynecologic surgery usage)



PM0902

Vaginal speculum with light source (Type A)



PM0903

Vaginal speculum with light source (Type B)



PM0904

Disposable anoscope with light source
(Type A)



PM0905

Disposable anoscope with light source
(Type B)



PM0906

Self-destroyed vaginal speculum



PM0907

Vaginal speculum medium screw type



PM0908

Screw fastener type



PM0909

Vaginal speculum side screw type



PM0910

Vaginal speculum fastener type



PM0911

Vaginal speculum rack type



PM0912

Vaginal speculum american type




PM0913

Vaginal speculum push-pull type



PM0914

Vaginal speculum butterfly type



PM0915

Vaginal speculum semilucent type



Surgical Blades

Gamma Radiation Sterilized
Barbon Steel material and Stainless Steel material to choose



Cod.	Specification	Material
Fitment No.3		
PM1001	10#	Barbon steel / Stainless steel
PM1002	11#	
PM1003	12#	
PM1004	12d#	
PM1005	13#	
PM1006	14#	
PM1007	15#	
PM1008	15d#	
PM1009	16#	
PM1010	17#	

Cod.	Specification	Material
Fitment No.4		
PM1011	18#	Barbon steel / Stainless steel
PM1012	19#	
PM1013	20#	
PM1014	21#	
PM1015	22#	
PM1016	23#	
PM1017	24#	
PM1018	24d#	
PM1019	25#	
PM1020	26#	
PM1021	27#	
PM1022	34#	
PM1023	36#	
PM1024	36d#	

Surgical Scissors



Cod.	Description	Length
PM1025	Surgical Scissors	12.5cm
PM1026		14cm
PM1027		16cm
PM1028		18cm
PM1029	Tissue Scissors	14cm
PM1030		16cm
PM1031		18cm
PM1032		20cm
PM1033	Suture Scissors	22cm
PM1034		14cm
PM1035	Bandage Scissors	14cm
PM1036		18cm

Cod.	Description	Length
PM1037	Surgical Hemostatic Forceps	12.5cm
PM1038		14cm
PM1039		16cm
PM1040		18cm
PM1041	Needle Holder	20cm
PM1042		22cm
PM1043		12.5cm
PM1044		14cm
PM1045		16cm
PM1046		18cm
PM1047		20cm
PM1048		22cm
PM1049		25cm

PM1050

Medical Forceps



PM1051

Cotton Balls Forceps




PM1052

Blood Forceps




PM1053

Dental Mirror




PM1054

Stainless steel surgical knife handle




PM1055

Disposable Scalpels



Stethoscope

PM1101- Single head stethoscope



PM1102- Single head stethoscope (Pediatric)



PM1103- Single head stethoscope with non-chill ring



PM1104- Single head stethoscope with non-chill ring (Pediatric)



PM1105- Sharp needle type stethoscope



PM1106- Sharp needle type stethoscope



PM1107- Dual head stethoscope



PM1108- Dual head stethoscope (Pediatric)



PM1109- Dual head stethoscope with non-chill ring




PM1110- Dual head stethoscope with non-chill ring (Pediatric)



PM1111- Neonatal and baby stethoscope



PM1112- Sprague rappaport stethoscope




PM1113- Golden sprague rappaport stethoscope




PM1114- Single tube sprague rappaport stethoscope



PM1115- Deluxe single tube sprague rappaport stethoscope



PM1116- Deluxe sprague rappaport stethoscope



**PM1117- Sprague
rappaport stethoscope
with clock**




**PM1118- Sprague
rappaport stethoscope
with "A" tube**



**PM1119- Deluxe dual
head stethoscope with
non-chill ring**



**PM1120- Dual head
stethoscope with inner-
spring binaural**



**PM1121- Deluxe dual
head stethoscope with
non-chill ring**



**PM1122- Special
single head
stethoscope**



**PM1123- Deluxe single
head stethoscope**



**PM1124- Cross-shaped
membrane lock single
head stethoscope**



**PM1125- Hump single
head stethoscope**




**PM1126- A/B Single
head stethoscope with
lock**



**PM1127- A/B Cardiology
stainless steel
stethoscope**



**PM1128- Classic
stainless steel
stethoscope**



**PM1129- Classic
stainless steel
stethoscope (Pediatric)**



**PM1130- Classic
stainless steel
stethoscope (Infant)**



**PM1131- Stainless steel
stethoscope (Neonate)**



**PM1132- Double diaphragm
cardiology stainless steel
stethoscope**



PM1133- Stainless steel single head stethoscope



PM1134- Cross-shaped membrane lock stainless steel stethoscope



PM1135- Three parts stainless steel stethoscope



PM1136- Bowl type stethoscope



PM1137- Flat head stethoscope



PM1138- Flat dual-purpose head stethoscope




PM1139- Deluxe single head with gold-plated stethoscope



PM1140- Deluxe single head stethoscope with "A" tubing




PM1141- Cardiology stainless steel stethoscope with gold-plated



PM1142- Classica stainless steel stethoscope with gold-plated



PM1143- Gold-plated cross-shaped membrane lock stainless steel stethoscope



PM1144- Cardiology single head stethoscope



PM1145- Cardiology single head stethoscope with gold-plated



PM1146- Matte/Bright classic single head stethoscope



PM1147- Acrylic stethoscope



PM1148- Teaching stethoscope



Sphygmomanometer


PM1149- Standard Aneroid Sphygmomanometer

- Zinc alloy gauge;
- Nylon standard cuff :54*14.5cm
- 2-Tube adult size latex bladder: 22x12cm
- Measure scope: 0-300mmHg;
- Accuracy: +/-3 mmHg, Sub-division: 2mmHg.



PM1150- Classic Aneroid sphygmomanometer

- Zinc alloy gauge;
- Nylon standard cuff with indicates print: 54x14.5cm
- 2-Tube adult size latex bladder: 22x12cm
- Measure scope: 0-300mmHg;
- Accuracy: +/-3 mmHg, Sub-division: 2mmHg.



PM1151- Economic Aneroid Sphygmomanometer

- Zinc alloy gauge;
- Nylon standard cuff: 54x14.5cm
- 2-Tube adult size PVC bladder: 22x12cm
- Measure scope: 0-300mmHg
- Accuracy: +/-3 mmHg, Sub-division: 2mmHg.



PM1152- Deluxe Aneroid Sphygmomanometer

- Zinc alloy gauge with brand ring, heavy weigh:145G
- Cotton standard cuff with unsticking
- 2-Tube adult size green latex bladder: 22x12cm
- Measure scope: 0-300mmHg;
- Accuracy: +/-3 mmHg, Sub-division: 2mmHg.



PM1153- Standard Kit Aneroid Sphygmomanometer With Stethoscope

- Zinc alloy gauge
- Cotton standard cuff with metal ring: 54x14.5cm
- 2-Tube adult size latex bladder: 22x12cm
- Measure scope: 0-300mmHg, Accuracy: +/-3 mmHg, Sub-division: 2mmHg
- Within single head stethoscope



PM1154- Digital Sphygmomanometer

- 60 group memories clock display
- Date and time display



PM1155- Digital Sphygmomanometer

- 60 group memories




PM1156- Digital Sphygmomanometer

- 60 group memories clock display



PM1157- Digital Sphygmomanometer

- 2x60 group memories date and clock display



Plaster




Cod.	Description	Size	Material
PM1201	Plastic Plaster	4x10cm/6x10cm/8x10cm/customizable	Plastic
PM1202	Fabric Plaster	4x10cm/6x10cm/8x10cm/customizable	Fabric
PM1203	Non-woven Plaster	4x10cm/6x10cm/8x10cm/customizable	Non-woven

Bandage




Cod.	Size	Material
PM1204	spot22mm/40x10mm/72x19mm/76x19mm/customizable	PE/PVC/EVA/PU
PM1205	spot22mm/40x10mm/72x19mm/76x19mm/customizable	Fabric/Elastic Fabric/Cotton
PM1206	56x19mm/72x19mm/76x19mm/customizable	PE/PVC/EVA/PU/Transparent Plastic

Cod.	Size	Material
PM1207	56X19mm	Cartoon PE/PVC
PM1208	56x19mm/72x19mm/76x19mm/customizable	Fabric/Elastic Fabric/Cotton



Cod.	Description	Size	Material
PM1209	Wound Dressing Plaster	6x7cm	Non-woven/PU
PM1210		7x9cm	
PM1211		7x10cm	
PM1212		6x12cm	
PM1213		9x10cm	

Cod.	Description	Size	Material
PM1214	Wound Dressing Plaster	9x15cm	Non-woven/PU
PM1215		9x20cm	
PM1216		9x25cm	
PM1217		9x30cm	
PM1218		9x35cm	




Cod.	Size	Color	Material	Packing Info.
Non-Woven Tape				
PM1219-1	1.25cm x 5m	Skin / White	Non-Woven	24rolls/box,30boxes/ctn
PM1219-2	2.5cm x 5m			12rolls/box,30boxes/ctn
PM1219-3	5cm x 5m			6rolls/box,30boxes/ctn
PM1219-4	7.5cm x 5m			6rolls/box,30boxes/ctn
PM1219-5	10cm x 5m			6rolls/box,30boxes/ctn
PE Tape				
PM1220-1	1.25cm x 5m	Transparent	PE	24rolls/box,30boxes/ctn
PM1220-2	2.5cm x 5m			12rolls/box,30boxes/ctn
PM1220-3	5cm x 5m			6rolls/box,30boxes/ctn
PM1220-4	7.5cm x 5m			6rolls/box,30boxes/ctn
PM1220-5	10cm x 5m			6rolls/box,30boxes/ctn
Silk Tape, Sawteeth				
PM1221-1	1.25cm x 5m	White	Silk	24rolls/box,30boxes/ctn
PM1221-2	2.5cm x 5m			12rolls/box,30boxes/ctn
PM1221-3	5cm x 5m			6rolls/box,30boxes/ctn
PM1221-4	7.5cm x 5m			6rolls/box,30boxes/ctn
PM1221-5	10cm x 5m			6rolls/box,30boxes/ctn
Zinc Oxide Plaster				
PM1222-1	1.25cm x 5m	Skin / White	Cotton	24rolls/box,30boxes/ctn
PM1222-2	2.5cm x 5m			12rolls/box,30boxes/ctn
PM1222-3	5cm x 5m			6rolls/box,30boxes/ctn
PM1222-4	7.5cm x 5m			6rolls/box,30boxes/ctn
PM1222-5	10cm x 5m			6rolls/box,30boxes/ctn
Zinc Oxide Adhesive Tape				
PM1223-1	1.25cm x 5m	Skin / White	Cotton	24rolls/box,30boxes/ctn
PM1223-2	2.5cm x 5m			12rolls/box,30boxes/ctn
PM1223-3	5cm x 5m			6rolls/box,30boxes/ctn
PM1223-4	7.5cm x 5m			6rolls/box,30boxes/ctn
PM1223-5	10cm x 5m			6rolls/box,30boxes/ctn
Zinc Oxide Tinplate Plaster				
PM1224-1	1.25cm x 5m	Skin / White	Cotton	12rolls/box,30boxes/ctn
PM1224-2	2.5cm x 5m			12rolls/box,30boxes/ctn
PM1224-3	5cm x 5m			6rolls/box,30boxes/ctn
PM1224-4	7.5cm x 5m			6rolls/box,30boxes/ctn

Cod.	Size	Color	Material	Packing Info.
Zinc Oxide Adhesive Tape				
PM1225-1	1/2" x 4 rolls	Skin / White	Skin / White	12IN. x 10yds
PM1225-2	1" x 3 rolls			
PM1225-3	2" x 2 rolls			
PM1225-4	3" x 1 rolls			
Zinc Oxide Plaster with Plastic Cover				
PM1226-1	1.25cm x 5m	Skin / White	Skin / White	24rolls/box,30boxes/ctn
PM1226-2	2.5cm x 5m			12rolls/box,30boxes/ctn
PM1226-3	5cm x 5m			6rolls/box,30boxes/ctn
PM1226-4	7.5cm x 5m			6rolls/box,30boxes/ctn
PM1226-5	10cm x 5m			6rolls/box,30boxes/ctn
Zinc Oxide Perforated Plaster				
PM1227-1	10cm x 5m	Skin	Cotton	1roll/box,80boxes/ctn
PM1227-2	10cm x 5m			10rolls/box,8boxes/ctn
PM1227-3	18cm x 5m			1roll/box,80boxes/ctn
PM1227-4	18cm x 5m			10rolls/box,8boxes/ctn
Kinesiology Tape				
PM1228-1	2.5cm x 5m	Skin/Pink/Blue/ others	95%Cotton+ 5%Spandex	6rolls/box,30boxes/ctn
PM1228-2	3.8cm x 5m			
PM1228-3	5cm x 5m			24rolls/box,30boxes/ctn
PM1228-4	7.5cm x 5m			12rolls/box,30boxes/ctn
PM1228-5	10cm x 5m			6rolls/box,30boxes/ctn
Sport Tape				
PM1229-1	2.5cm x 5m	White	Cotton	12rolls/box,30boxes/ctn
PM1229-2	5cm x 5m			6rolls/box,30boxes/ctn
PM1229-3	3.8cm x 5m			6rolls/box,30boxes/ctn
Non-Woven Cohesive Bandage				
PM1230-1	2.5cm x 4.5m	Printing	Non-Woven	24rolls/box,24boxes/ctn
PM1230-2	5cm x 4.5m			12rolls/box,24boxes/ctn
PM1230-3	7.5cm x 4.5m			12rolls/box,16boxes/ctn
PM1230-4	10cm x 4.5m			12rolls/box,8boxes/ctn
Cotton Cohesive Bandage				
PM1231-1	2.5cm x 4.5m	Red/Black/ Green/White/ Others	Cotton	24rolls/box,24boxes/ctn
PM1231-2	5cm x 4.5m			12rolls/box,24boxes/ctn
PM1231-3	7.5cm x 4.5m			12rolls/box,16boxes/ctn
PM1231-4	10cm x 4.5m			12rolls/box,8boxes/ctn



Cod.	Model	Product Dimension (mm)	Strip Size (mm)	Material
PM1301	K3	90 x 12	60 x 3.0	ABS
PM1302	M8SX	70 x 20	60 x 4.0	ABS
PM1303	LK4-1	70 x 20	60 x 4.0	ABS
PM1304	WE-H	70 x 20	60 x 4.0	ABS
PM1305	WD	70 x 22	60 x 3.5	ABS

Cod.	Model	Product Dimension (mm)	Strip Size (mm)	Material
PM1306	LK2	75 x 20	60 x 4.0	ABS
PM1307	CARE	70 x 35	60 x 4.0	ABS
PM1308	AP	90 x 16	80 x 4.0	ABS
PM1309	RY	90 x 18	80 x 4.0	ABS




Cod.	Model	Product Dimension (mm)	Strip Size (mm)	Material
PM1310	Multi-3	70 x 40	60 x 3.0	ABS
PM1311	Multi-3	70 x 44	60 x 4.0	ABS
PM1312	TZYK Multi-3	70 x 45	60 x 4.0	ABS
PM1313	ZX Multi-3	70 x 45	60 x 3.0	ABS
PM1314	Multi-2	93 x 27	4.0	ABS

Cod.	Model	Product Dimension (mm)	Strip Size (mm)	Material
PM1315	Multi-3	93 x 27	4.0	ABS
PM1316	Multi-4	93 x 45.5	4.0	ABS
PM1317	Multi-5	93 x 45.5	4.0	ABS
PM1318	Multi-7	93 x 45.5	4.0	ABS
PM1319	Multi-8	93 x 45.5	4.0	ABS

Disposable Drainage Bag Series

- The Disposable Drainage Bag Series is used to collect and measure the volume of urine from a patient, either during the course of normal urination or while the patient is catheterized.
- The length of tube can be customized according to customer requirements.



Cod.	Description	Volume	Specification
PM1401	Urine Bag without outlet	1500ml / 2000ml	---
PM1402	Urine Bag with Push-Pull Valve	1000ml / 1500ml / 2000ml	Anti-Reflux Flutter Valve, Push-Pull Valve
PM1403	Urine Bag with Cross Valve	1000ml / 1500ml / 2000ml	Anti-Reflux Flutter Valve, Cross Valve
PM1404	Urine Bag with Twist Turn Valve	1000ml / 1500ml / 2000ml	Twist Turn Valve
PM1405	Urine Bag for Leg	500ml / 750ml	Anti-Reflux Flutter Valve, Push-Pull Valve, Two Pre-assembled Straps
PM1406	Luxury Urine Bag I	2000ml	Rope Hanger, Adapter with Cap, Anti-Reflux Valve, Sampling Port
PM1407	Luxury Urine Bag II	2000ml	Rope Hanger, Adapter with Cap, Anti-Reflux Valve
PM1408	Luxury Urine Bag III	2000ml	Rope Hanger, Adapter with Cap, Anti-Reflux Valve, Outlet Tube with Slide Clamp
PM1409	Padiatric Urine Bag	---	---

PM1501 Sterile Gauze Swab

- 100% cotton, high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 11,12,13,15,17, 20, 22 threads etc.
- With or without x-ray detectable threads
- Folded edge or unfolded edge
- 5x5cm,7.5x7.5cm,10x10cm,10x20cm, 20x20cm etc.
- 8ply,12ply,16ply, 24ply, 32ply etc.
- Sterile: gamma, EO, steam
- Pouch: paper+paper, paper+film
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1503 Cutting Gauze

- 100% cotton, high absorbency and softness
- Cotton yarn of 21's, 32's
- Mesh of 11, 12, 13, 15, 17, 20, 22 threads etc.
- 10x10cm,10x20cm, 20x20cm, 20x30cm, 30x30cm etc.
- For beauty,cosmetics,dressing wound, scrubbing of medical equipment and apparatus
- Can be customized according to customers' sample, drawing, customer requirements and any specifications




PM1505 Non-Sterile Lap Sponge / Abdominal Pad

- 100% cotton, high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 11, 12, 13, 15, 17, 20, 22 threads etc.
- With or without x-ray detectable threads
- Pre -washed or non-washed
- 20x20cm, 22.5x22.5cm, 30x30cm, 40x40cm, 45x45cm, 50x50cm etc.
- 4ply, 6ply, 8ply, 12ply, 16ply, 20ply etc.
- Ribbon: blue, green, etc.
- Can be customized according to customers' sample, drawing, customer requirements and any specifications




PM1502 Non-Sterile Gauze Swab

- 100% cotton, high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 11,12, 13,15,17, 20, 22 threads etc.
- With or without x-ray detectable threads
- Folded edge or unfolded edge
- 5x5cm, 7.5x7.5cm, 10x10cm, 10x20cm, 20x20cm etc.
- 8ply,12ply,16ply, 24ply, 32ply etc.
- 50pcs/pack ,100pcs/pack, 200pcs/pack etc..
- Can be customized according to customers' sample, drawing, customer requirements and any specifications




PM1504 Sterile Lap Sponge / Abdominal Pad

- 100% cotton, high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 11, 12, 13, 15, 17, 20, 22 threads etc.
- With or without x-ray detectable threads
- Pre-washed or non-washed
- 20x20cm, 22.5x22.5cm, 30x30cm, 40x40cm, 45x45cm, 50x50cm etc.
- 4ply, 6ply, 8ply, 12ply, 16ply, 20ply etc.
- Sterile: gamma, EO, steam
- Ribbon: blue, green, etc.
- Pouch: 1's, 2's, 5's/csk package blister packing
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1506 Sterile And Non-Sterile Gauze Ball

- 100% cotton,high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 17, 20 threads
- With or without x-ray detectable threads
- Pre-washed or non-washed
- 8x8cm, 9x9cm, 15x15cm, 18x18cm, 20x20cm, 25x30cm etc.
- Sterile: gamma, EO, steam
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1507 Gauze Roll

- 100% cotton high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 11, 12, 13, 15, 17, 20, 22 threads etc.
- With or without x-ray detectable threads
- 1 ply, 2ply, 4ply, 8ply
- 36"x5m, 36"x5yds, 36"x50m, 36"x50yds, 36"x100m, 36"x100yds
- Packing: 1roll /blue Kraft paper or poly-bag
- Can be customized according to customers' sample, drawing, customer requirements, and any specifications



PM1508 Gauze Bandage

- 100% cotton, high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 11, 12, 13, 15, 17, 20, 22 threads etc.
- Width: 5cm, 7.5cm, 14cm, 15cm, 20cm
- length: 3m, 3yds, 4m, 4yds, 5m, 5yds, 7m, 7yds, 10m, 10yds.
- 36"x5m, 36"x5yds, 36"x50m, 36"x50yds, 36"x100m, 36"x100yds
- Packing: 10roll /pack, 12roll /pack (non sterile)
1roll /pouch (sterile)
- Sterile: gamma, EO, steam
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1509 Zigzag Gauze / Pillow Gauze

- 100% cotton, high absorbency and softness
- Cotton yarn of 21's, 32's, 40's
- Mesh of 10, 13, 17, 20 threads etc.
- Width: 36" (90cm), 48" (120cm)
- Length: 5m, 10m, 20m, 30m, 50yards, 50m, 100yards, 100m, 200m, 500m etc.
- Layer: 1ply, 2ply, 4ply
- With or without x-ray detectable threads
- Packing: 1 roll/blue Kraft paper or PE bag
- Folded in roll, Zigzag, Pillow or cut
- Can be customized according to customers' sample, drawing, customer requirements, and any specifications



PM1510 Cotton Wool Roll/Cotton Ball

- 100% cotton, high absorbency and softness
- Weight of roll: 25g, 50g, 100g, 150g, 200g, 250g, 400g, 500g, 1000g, etc.
- Packing for roll: 1 roll/blue Kraft paper or PE bag
- Weight of ball: 0.15g, 0.25g, 0.4g, 0.5g, 0.8g, 1.0g, 2.0g, 3.0g, etc.
- packing for ball: 100g/bag, 250g/bag, 500g/bag
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1511 Face Mask

- Disposable surgical /medical face mask
- Layer: 1ply/2ply/3ply/4ply
- Adjustable Nose Clip
- With tie or with ear loop
- Outer layer: 20gsm/M2 PP Non Woven
- Middle layer: Melt Blow PP Fiche Paper BFE>95%/99% Available
- Inner layer: 25gsm/M2 PP Non Woven
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1512 Non Woven Swab

- Sterile or non sterile non woven swab
- Material: non woven
- Sterilization: EO, gamma
- Layer: 4ply, 6ply, 8ply, 12ply etc.
- Size: 5x5cm, 7.5x7.5cm, 10x10cm or customized
- Weight: 30g, 35g, 40g, 50g, etc
- X-ray: With or without x-ray
- Feature: high absorbency and softness
- Packing:
Sterile: 1's, 2's, 10's, etc packed in pouch
Non-sterile: 100pcs/pack, or as your required
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1513 Surgical/Operation Gown

- Disposable surgical /medical gown
- Material:
 1. PP/S (100% Polypropylene Spun bond Non woven Fabric)
 2. MS (Polypropylene Spun bond Non woven Fabric + Melt blown Non woven fabric + Polypropylene Spun bond Non woven Fabric)
 3. P/S+PE Film
- Size: S(115cm*137cm) M(120cm*140cm) L(130cm*150cm) XLL(140cm*160cm) or any other customized sizes
- Weight: 18-60gsm is available(as your request)
- Color: Blue, white, green, purple
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1515 Non Woven Bed Sheet

- Disposable surgical /medical n-w bed sheet
- Material: SBPP/SMS/PP+PE
- With elastic/ties
- Weight: 15-45gsm as your requirement
- Size: 170x230cm, 120x220cm, 100x180cm or as your requirement
- Color: white, blue, green, pink, etc.
- Can be customized according to customers' sample, drawing, customer requirements and any specifications




PM1514 Non Woven Cap

- Disposable surgical /medical n-w cap
- Material: PP or SMS
- Bouffant Cap/Doctor Cap/Clip Cap
- With Tie-on or Elastic Band
- Color: Blue, white, green, purple
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1516 Non Woven Bed Cover

- Disposable surgical /medical n-w bed cover
- Material: SBPP/SMS/PP+PE
- With elastic/ties
- Weight: 15-45gsm as your requirement
- Size: 80x210cm, 100x210cm, 110x240cm or as your requirement
- Color: white, blue, green, pink, etc.
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



PM1517 Non Woven Shoe Cover

- Disposable surgical /medical n-w shoe cover
- Material: PP or PE
- With/without anti-skid
- Weight: 25gsm, 30gsm, 35gsm
- Size: 15x40cm, 17x38cm, 17x41cm, 17x45cm
- Color: Blue, white, green, etc.
- Can be customized according to customers' sample, drawing, customer requirements and any specifications



Alcohol Product

Alcohol Pad(70% IPA)



Cod.	Spec.	Tissue Size
PM1601	5x5cm	6x3cm 2ply
PM1602	5x5cm	6x6cm 4ply
PM1603	6.2x7cm	9x4.2cm 2ply
PM1604	7.5x10cm	15x20cm 8ply

Povidone-Iodine Prep Pad



Cod.	Spec.	Tissue Size
PM1605	5x5cm	6x3cm 2ply
PM1606	5x5cm	6x6cm 4ply
PM1607	6.2x7cm	9x4.2cm 2ply
PM1608	7.5x10cm	15x20cm 8ply

Alcohol+Chg Pad



Cod.	Spec.	Tissue Size
PM1609	5x5cm	6x3cm 2ply
PM1610	5x5cm	6x6cm 4ply
PM1611	6.2x7cm	9x4.2cm 2ply
PM1612	7.5x10cm	15x20cm 8ply

BZK Pad



Cod.	Spec.	Tissue Size
PM1613	5x5cm	6x3cm 2ply

Clean Rubber Wiper



Cod.	Spec.	Tissue Size
PM1614	5x5cm	6x3cm 2ply

Alcohol Pad(70% IPA)



Cod.	Spec.	Tissue Size
PM1615	10x10cm	15x15cm 6ply

Saline Swab



Cod.	Spec.	Tissue Size
PM1616	5x5cm	6x3cm 2ply

Antiseptic Wipe



Cod.	Spec.	Tissue Size
PM1617	6x5cm	13x20cm 32ply

Soap Wipe



Cod.	Spec.	Tissue Size
PM1618	6x5cm	13x20cm 32ply

Sterilizing Cotton Swabs-(Alcohol)



Cod.	Spec.	Stick
PM1619	4"xØ10mm absorbent cotton	4", plastic

Sterilizing Cotton Swabs-(Povidone)




Cod.	Spec.	Stick
PM1620	4"xØ10mm absorbent cotton	4", plastic

Sterilizing Cotton Swabs-(Alcohol+Chg)




Cod.	Spec.	Stick
PM1621	4"xØ10mm absorbent cotton	4", plastic

PM1701 Colostomy Bag Type A




PM1702 Colostomy Bag Type B




PM1703 Umbilical Clamp



PM1704 Umbilical Clamp Cutter



PM1705 Sponge Brush Type A



PM1706 Sponge Brush Type B



PM1707 Oral cavity kit



PM1708 Dressing Kit Type A




PM1709 Dressing Kit Type B




PM1710 Dressing Kit Type C



PM1711 Tracheostomy care kit



PM1712 Pap smear kit Type A



PM1713 Pap smear kit Type B



PM1714 Hand wash brush




PM1715 Hand wash brush



PM1716 Shaver Type A



PM1717 Shaver Type B



PM1718 Digital Quartz Timer



PM1719 Mechanical Timer



PM1720 Digital Timer Type A



PM1721 Digital Timer Type B



PM1722

Funnels with Perforated Filter Dish
Glazed, the o.d. is from 20mm to 300mm



PM1723

Mortars with Spout, Pestles Glazed
The o.d. is from 60mm to 305mm



PM1724

Porcelain crucibles
high / medium / low wall with lid



PM1725

Porcelain Basins, round bottom with spout
The o.d. is from 47mm to 118mm




PM1726

Stir Bar Egg shaped



PM1727

Stir Bar Octagon shaped



PM1728

Stir Bar
Octagon shaped with a center molded pivot ring







Cod.	Description	Specification
PM1729	Wooden Tongue Depressor 150x18x1.6mm	bulk packing
PM1730		individual packing
PM1731		for pediatric



Cod.	Description	Specification
PM1732	Triangles	50mm
PM1733		63mm
PM1734		75mm



Cod.	Description	Specification
PM1735	Ceramic wire gauze	10x10cm
PM1736		12.5x12.5cm
PM1737		15x15cm
PM1738	Ceramic wire gauze	18x18cm
PM1739		20x20cm
PM1740		25x25cm




Cod.	Description	Specification
PM1741	Identification bracelet for adult	write-on
PM1742		with card insert
PM1743	Identification bracelet for kids	write-on
PM1744		with card insert



Cod.	Description	Specification
PM1745	Serological pipette sterile container	Stainless steel material
PM1746	Petri dish sterile container	Stainless steel material
PM1747	Medical Paper	50mm x 30m
PM1748		60mm x 30m
PM1749		50mm x100mm-300p
PM1750	Laboratory Thermometer	110mm x 140mm-200p
PM1751		Red spirit-filled
PM1752		Mercury-filled

Cod.	Description	Specification
PM1753	Ear Syringe TPE Material	30ml
PM1754		60ml
PM1755		90ml
PM1756	Ear Syringe PVC Material	30ml
PM1757		60ml
PM1758		90ml
PM1759	Pipette Filter	Red / Green / Black
PM1760	Silicon Cap for Reagent Bottle	Ø12(9-14)
PM1761		Ø10(8-12)



Cod.	Description	Capacity	Specification
PM1801	Phlebotomy Sharps Containers	1qt	Slide Lid
PM1802		1qt	Hinged Lid
PM1803		2qt	Slide Lid
PM1804		2qt	Rotor-opening Lid

Accessories:



Holders



Holders



Brackets



Cod.	Description	Capacity	Specification
PM1805	In Room Sharps Containers Type A	5qt	Mailbox-Lid
PM1806		2gal	Mailbox-Lid
PM1807		3gal	Mailbox-Lid
PM1808	In Room Sharps Containers Type B	5.4qt	Mailbox-Lid
PM1809		2gal	Mailbox-Lid
PM1810		3gal	Mailbox-Lid

Accessories:



Locking Wall Bracket



Locking Wall Cabinet



Glove Box Holder for Wall Cabinets



Cod.	Description	Capacity	Specification
PM1811	Multi-Purpose Sharps Containers	1gal	Slide Lid
PM1812		2gal	Slide Lid
PM1813		3gal	Slide Lid

Cod.	Description	Capacity	Specification
PM1814	Multi-Purpose Sharps Containers	6gal	Slide Lid
PM1815		8gal	Slide Lid



Cod.	Description	Capacity	Specification
PM1816	Type A	5qt	Rotor-opening Lid
PM1817		2gal	Rotor-opening Lid
PM1818		3gal	Rotor-opening Lid
PM1819	Type B	2gal	Hinged&Rotor-opening Lid
PM1820		3gal	Hinged&Rotor-opening Lid

Cod.	Description	Capacity	Specification
PM1821	Type C	8gal	Hinged&Rotor-opening Lid
PM1822	Type D	8qt(2gal)	Vertical-Drop Lid
PM1823		14qt	Vertical-Drop Lid
PM1824		6gal	Vertical-Drop Lid

Transportable Containers

- Simple, one-handed operation opens and closed the one quart disposal opening
- Leak-resistant screw cap offers secure closure between uses when properly closed
- Large disposal opening easily accommodates small syringes and other sharps waste



Cod.	Description	Capacity	Specification
PM1825	Type A	0.2L	Slide Lid
PM1826	Type B	0.3L	Hinged Lid
PM1827	Type B	0.3L	Hinged Lid

Cod.	Description	Capacity	Specification
PM1828	Type C	1QT	Hinged Lid
PM1829	Type C	1.5QT	Hinged Lid

Sharps Containers

- Made of polypropylene and in square shape
- Leak proof and puncture resistant, confirm to EN ISO 23907 compliance
- Slide Lids or Rotating lids for the disposal of sharp and cutting waste in different sizes
- Operates temporary closure and final closure
- Ergonomic handle for the more spacious models
- Maximum filling level indicator
- Fitted with a universal trolley hanging system



Cod.	Description	Capacity	Specification
PM1830	Sharps Containers	2L	Slide Lid
PM1831		4L	Slide Lid
PM1832		5L	Slide Lid

Cod.	Description	Capacity	Specification
PM1833	Sharps Containers	15L	Slide Lid
PM1834		23L	Slide Lid



Cod.	Description	Capacity	Specification
PM1835	Sharps Containers	3L	Rotor-opening Lid, Handle
PM1836		5L	Rotor-opening Lid, Handle
PM1837		6L	Rotor-opening Lid, Handle

Cod.	Description	Capacity	Specification
PM1838	Sharps Containers	10L	Mailbox-Lid

Transportable Containers

- Made of polypropylene and in round shape
- Mechanism for unhooking insulin needles
- Leak proof and puncture resistant, confirm to EN ISO 23907 compliance
- Operates Temporary closure and final closure
- Maximum filling level indicator
- Various supports available



Cod.	Description	Capacity	Specification
PM1839	Sharps Containers	0.8L	Rotor-opening Lid
PM1840		2.8L	Rotor-opening Lid
PM1841		3.6L	Rotor-opening Lid


Cod.	Description	Capacity	Specification
PM1842	Sharps Containers	6.0L	Rotor-opening Lid, Handle
PM1843		6.0L	Slide Lid, Handle



Cod.	Description	Capacity	Specification
PM1844	Sharps Containers	0.7L	Hinged Lid
PM1845		1.0L	Hinged Lid
PM1846		1.5L	Hinged Lid
PM1847		2.0L	Hinged Lid
PM1848		3.0L	Hinged Lid
PM1849		6.5L	Hinged Lid
PM1850		8.0L	Hinged Lid

Needle Counter Box

- Medical grade high impact plastic to prevent accidental needle sticks
- High-density foam holds needles in place securely detachable halves allow the same needle counter to be used in two separate areas
- Snap closure keeps needle counter contents in one concise package
- Adhesive tabs to hold each half securely to any surface and lock the halves closed for safe disposal Latex free




Cod.	Description	Capacity	Specification
PM1851	Needles Counter Box	20-40L Count Double Magnet	11.5x5.5x1.5cm
PM1852		20L count Foam Block/Magnet	11.5x6.5x3.3cm

Cod.	Description	Capacity	Specification
PM1853	Needles Counter Box	30L Count Foam Strip/Adhesive	12.2x10.5x1.7cm
PM1854		40L Count Foam Block/Adhesive	12.2x10.5x3.5cm

Disposable Surgical Pack




Cod.	Description	Capacity	Specification
PM1844	Basin	250cc	Cup, Denture W/O Lid
PM1845		500cc	Bowl, 16 OZ
PM1846		1000cc	/
PM1847		2500cc	/
PM1848		6000cc	Ring
PM1849		6500cc	Ring
PM1846	Lid	250cc	For Denture Cup
PM1847		500cc	For 16OZ Bowl
PM1848		1000cc	For 1000cc Basin
PM1850		2500cc	For 2500cc Basin



Cod.	Description	Capacity	Specification
PM1844	Tray	/	Quarter Shallow 4.7x5.4x1"
PM1845		/	Quarter Deep 4.7x5.4x2"
PM1846		/	Half Shallow
PM1847		/	Half Deep
PM1848		/	Half, 2 Comp 4.6x9.6x1.8"
PM1849		/	Full, 3 Comp 9.6x10x2" Plastic
PM1846		/	Full Deep



Cod.	Description	Capacity	Specification
PM1844	Guidewire Basin	2000cc	Guidewire Basin
PM1845	Placenta Basin	2500cc	Bowl, Guidewire No Lid Ridge Bulk
PM1846		3000cc	
PM1847	Lid	3000cc	For Placenta Basin
PM1848	Basin	8000cc	Rect 10x13-3/4x5" H
PM1849	Emesis Basin	700cc	Emesis Basin, Blue 24OZ
PM1846		700cc	Orange Emesis



Cod.	Description	Capacity	Specification
PM1844	Tray	/	Small Base
PM1845		/	Large Deep
PM1846		/	Large Mayo
PM1847		/	Large Base
PM1848		/	Small Base
PM1846		/	Large Mayo



宁波普瑞供应链服务有限公司

PRIMARY
SCS

NINGBO PRIMARY SUPPLY CHAIN SOLUTION CO.,LTD.

6F, BLOCK C, 89 CUIBAI ROAD, NINGBO 315010, CHINA

Tel: 86-574-8386 7806

Cell: 86-139 6782 7576

Fax: 86-574-8386 7808

E-mail: tracy@primaryscs.com

tracy399@163.com